

OSSTEM[®]
IMPLANT

2020 – 21
**KIT PRODUCT
CATALOG**

CE 2460
PG IM22
ENG_ver.

OSSTEM[®]
IMPLANT

Osstem Implant
2020-21 Comprehensive Catalog

Planning/Editing Design Center Brand Design Team
Supervision PM Headquarters, Implant Lab
Production/Distribution PM Headquarters
Date of Publication 02.2020
Publisher OSSTEM IMPLANT Co., Ltd. 3, Magokjungang 12-ro,
Gangseo-gu, Seoul 07789, (Republic of) Korea
Phone +82.2.2016.7000
Fax +82.2.2016.7001
www.osstem.com

KIT PRODUCT CATALOG

003	INTRODUCTION
012	CONTENTS
022	OSSTEM KIT
210	User Manual

"Cutting edge
technology and
superior quality"

Making products that dentists are able to
trust and are satisfied with:
This is our mission at **OSSTEM IMPLANT**

We are forever grateful
to all of our
customers for their
unwavering support to
OSSTEM IMPLANT

Osstem, South Korea's first implant manufacturer, has achieved steady growth thanks to the support and love from its customers. Osstem has put a lot of effort into continuous investment in R&D and quality innovation in order to provide products that customers are looking for and satisfied with. Based on this, it has become the No. 1 implant company in Asia Pacific region and No. 4 in the world. Moreover, it was ranked No. 1 for global fixture sales from 2017 to 2019 for 3 consecutive years and became the global provider of the implants most used by the customers all around the world.

In this 2020-21 product catalog, you can see a variety of products at a glance, including not only the implant products of Osstem's differentiated technology but also the digital dentistry products such as Oneguide the implant surgical guide, scanners, milling machines, CAD/CAM, etc. We have invested numerous efforts and time in the configuration and design of this catalog so that customers do not have any inconvenience of finding and ordering the products they need. The fixtures and abutments are listed to make it easy to understand the diameter, length, and functional behavior, in sequence that customers make a judgement for purchase.


The product type and code are displayed to help with accurate ordering. We have added the product thumbnail pages to view the components at a glance and detailed information pages to describe the functions of each component for enhanced user understanding. For GBR products, shape, size, and capacity of each product are described in detail for easy ordering as well. In addition, the release date and time are indicated for all products so that customers can easily distinguish new products from existing products for purchase. In terms of design, we applied high-quality product images to aid ordering without looking at the actual product, and improved user convenience by applying representative colors to facilitate classification by product category.

We hope that this 2020-21 product catalog will help you effectively find and purchase all the products you need for your dental practice. Osstem Implant will continue to strive to create greater customer value as a partner to help dentists provide better care. Thank you.

CEO of OSSTEM IMPLANT
Tae-Kwan Eom



Worldwide & History



Headquarters
Overseas Subsidiary
Distributor

EMEA

GERMANY
RUSSIA
UKRAINE
FRANCE
HUNGARY
SWITZERLAND
NORWAY
FINLAND
UK
POLAND
ROMANIA
CZECH
SLOVAKIA
SERVIA
PORTUGAL
SPAIN
ITALY
CROATIA
GREECE
LATVIA
ESTONIA
LEBANON
TUNISI
MACEDONIA
SLOVENIA
KOSOVO
BULGARIA
GEORGIA
EGYPT
SOUTH AFRICA
ALBANIA
TAJIKISTAN

ASIA / OCEANIA

KOREA
JAPAN
CHINA
CHINESE TAIPEI
VIETNAM
BANGLADESH
HONG KONG, CHINA
SAUDI ARABIA
AUSTRALIA
NEW ZEALAND
TURKEY
PAKISTAN
KUWAIT
U.A.E
IRAN
CAMBODIA
PAPUA NEW GUINEA
MYANMAR
SRI LANKA
OMAN
JORDAN
PALESTINE
UZBEKISTAN

N/S.AMERICA

CANADA
USA
MEXICO
CHILE
COSTA RICA
PERU
BRAZIL

1997

- 01 Established Osstem R&D System
- 12 Launched "Doobun" (Health insurance claiming software)

2000

- 06 Developed and launched "Hidaro" (total dental clinic management software)
- 12 Acquired Sumcomprehensive Dental Materials (South Korea's first implant manufacturer)

2001

- 01 Obtained CE-0434 certification
- 03 Established AIC Training Center

2002

- 01 Established Osstem Implant Research Center
- 08 Obtained US FDA certification

2003

- 07 Established the Information System Research Institute

2006

- 03 Changed company name to Osstem Implant Co., Ltd.
- 09 Established a subsidiary in the U.S. (HIOSSSEN), and set up the manufacturing facility
- 12 Completed the first-phase establishment of overseas subsidiaries (12 countries)

2007

- 02 Listed on KOSDAQ and began trading
- 11 Won the "10 Million Dollar Export Tower" on Trade Day

2008

- 01 Established Osstem Bone Science Research Institute
- 07 Won the Grand Prize of the 2008 Korea Health Industry Awards by the Ministry of Health, Welfare and Family Affairs

2010

- 03 Launched TSIII SA implant
- 06 Launched TSIII HA implant

2011

- 06 Osstem Implant Research Institute selected as an Advanced Technology Center (ATC) by the Ministry of Trade, Industry and Energy
- 07 Selected as 2011 World Champ company by KOTRA
- 12 Selected as Current World-Class Product by the Ministry of Knowledge Economy

2012

- 06 Launched TSIII CA implant
- 07 Established the Medical Equipment Research Institute

2013

- 01 Launched xenograft "A-Oss"
- 09 Launched "K3 unit chair"

2014

- 05 Launched impression material "Hysil"
- 08 Launched whitening material "BeauTis"

2015

- 03 Established Osstem Pharma Co., Ltd.
- 12 Awarded the "50 Million Dollar Export Tower" on Trade Day

2016

- 01 Established VUSSEN Co., Ltd.
- 02 Released TSIII BA
- 03 Acquired Cardiotec Co., Ltd.
- 04 Launched the dental clinic interior design business
- 06 Released TSIII SOI
- 08 Acquired Hubit Co., Ltd.
- 11 Launched "OneGuide"

2017

- 12 Won the Presidential Award at 2017 Government Commendation for Job Creation
- 11 Won the '2018 SW Enterprise Quality Award' by Ministry of Science and Technology
- 12 Won the "100 Million Dollar Export Tower" on Trade Day


2018

2019

- 08 Opened manufacturing corporation in Yancheng, China
- 10 Established a subsidiary in Brazil (23 subsidiaries in 26 countries in operation)
- 12 Awarded the Brand Top, Industrial Service, Presidential Citation, Prime Minister Citation, and KITA Citation on 56th Trade Day

2020

- 01 Launched "OneClick" the electronic chart for dental clinics
- 02 Established "DenAll", the comprehensive dental portal
- 07 Headquarters relocated to Magok, Seoul
- 08 No.1 seller of fixture for 3 consecutive years (2017~2019)



Next-generation submerged type implant with an Internal hex 15° tapered connection structure

- Connection - **Regular only** (2.1hex single platform)
 - Strength intensified due to a narrower and deeper connection
 - Reduced prosthetic errors and inventory burden with no variation of the product (Mini/Regular)
- Abutment holding system applied to enable screw fastening with one hand
- Excellent initial stability in soft bone with smaller threads in the upper section
- Corkscrew thread & cutting edge
 - Superior self-threading effect for easy placement path adjustment
 - Enhanced initial stability in soft bone and consistent placement torque according to the drill diameter
- Available surface types - BA

Submerged type implant with an internal hex 11° tapered connection structure

- Connection - **Mini / Regular**
- Excellent initial stability in soft bone with smaller threads in the upper section
- Corkscrew thread & cutting edge
 - Superior self-threading effect for easy placement path adjustment
 - Enhanced initial stability in soft bone and consistent placement torque according to the drill diameter
- Various body shape options available to match the patient's bone quality and clinical condition
 - TSII (straight body) : Easy to adjust placement depth
 - TSIII (1.5° tapered body) : Excellent initial stability needed for immediate loading even in soft bone
 - TSIV (6° tapered body) : Specifically designed for use in maxillary sinus and soft bone, providing excellent initial stability
- Available surface types - SA / CA / BA / SOI

Non-submerged type implant with an internal octa 8° tapered connection based on 1st stage surgery

- Connection - **Regular / Wide**
- Corkscrew thread & cutting edge
 - Superior self-threading effect for easy placement path adjustment
 - Enhanced initial stability in soft bone and consistent placement torque according to the drill diameter
- Various body shape options available to match the patient's bone quality and clinical condition
 - SSII (straight body) : Easy to adjust placement depth
 - SSIII (1.5° tapered body) : Excellent initial stability needed for immediate loading even in soft bone
- Available surface types - SA / CA / BA


Submerged type implant with an external hex connection structure

- Connection - **Mini / Regular / Wide / Wide PS**
- Corkscrew thread & cutting edge
 - Superior self-threading effect for easy placement path adjustment
 - Enhanced initial stability in soft bone and consistent placement torque according to the drill diameter
- Various body shape options available to match the patient's bone quality and clinical condition
 - USII (straight body) : Easy to adjust placement depth
 - USIII (1.5° taper body) : Excellent initial stability needed for immediate loading even in soft bone
 - USIV (6° taper body) : Specifically designed for use in maxillary sinus and soft bone, providing excellent initial stability
- Available surface types - SA / CA / BA / SOI

OSSTEM[®] Implant

The key factor in providing implant treatment safely and efficiently is surface technology.
OSSTEM IMPLANT proudly presents its cutting-edge surface technology

Surface feature



Optimized Surface through Acid Treatment

- Ra 2.0~3.0 \square surface roughness
(Note : The roughness in the upper 0.5mm part is Ra 0.5~0.6 \square)
- Consistent surface micro-pits of 1~3 \square
- Surface area increased by 46% compared to RBM treated implants

In-vitro and In-vivo Bone Response

- Osteoblast differentiation and ossification improved by 20% compared to RBM treated implants
- Initial bone reaction performance in big animal model (mini-pig)
 - Initial stability (RT, 4 weeks) improved by 48% compared to RBM treated implants
 - Ossification (BIC, 4 weeks) improved by 20% compared to RBM treated implants

Super-hydrophilic SA surface immersed in a calcium solution

- Same surface morphology as SA surfaces
- Surface reaction activated by immersing in a calcium solution (CaCl₂)
- Increased new bone formation area with excellent blood wettability
- Bone response improved in early osseointegration stage compared to standard SA surface

In-vitro and In-vivo Bone Response

- Protein and cellular adhesion tripled compared to SA surfaces
- Initial cellular differentiation (7 days) improved by 19% compared to SA surfaces
- Initial stability (RT, 4 weeks) improved by 34% compared to SA surfaces
- Ossification (BIC, 4 weeks) improved by 26% compared to SA surfaces

Premium low crystalline nano-HA coated SA surface

- 10nm ultra-thin HA coating
- SA surface (Ra 2.0~3.011 \square) coated with HA
- Dual functions of titanium and HA
 - HA is naturally resorbed during ossification

In-vitro and In-vivo Bone Response

- Advantages of both SA and HA surfaces
 - SA's ability to maintain an optimal surface
 - HA's ability to form high quality initial bone even in bone of poor quality
- Ossification (BIC) improved by 26% compared to SA surfaces
- Applicable to all types of bone quality

Next-generation surface with hemostatic effect and pH control feature

- Activation of blood clot formation
- Prevention of carbon adsorption in air
- Same surface roughness (Ra 2.0~3.0 \square) as SA surfaces
- Superior blood wettability with super hydrophilic surface

In-vitro and In-vivo Bone Response

- Protein and cellular adhesion increased by 130 times compared to SA surfaces
- Initial stability (RT, 4 weeks) improved by 57% compared to SA surfaces
- Surface with the shortest duration of treatment







































































KIT Contents 1/4

022 One Guide KIT	023 Tissue Punch	023 Flattening Drill	024 Initial Drill	024 Initial Drill (Short type)	055 MS Narrow Ridge Fixture Driver	056 Fixture Driver (TS)	056 Fixture Driver (TS, Stopper Type)	056 Adapter	057 Driver Separator
024 OneGuide Twist Drill (Ø2.2)	025 OneGuide Drill	026 OneGuide Taper Cortical Drill	026 NoMount Driver	027 Fixture Driver	057 OneMS Path Drill	057 OneMS Lance Drill	060 One CAS KIT	061 OneCAS Twist Drill (Ø2.2)	061 OneCAS Drill
027 Fixture Driver (Stopper type)	028 OneGuide SS NoMount Driver	028 OneGuide SS Fixture Driver	029 OneGuide US NoMount Driver	029 OneGuide US Fixture Driver	062 OneCAS Stopper	062 Depth Gauge	063 Hydraulic Membrane Lifter	063 Bone Carrier Head	063 Bone Carrier
030 OneGuide KS NoMount Driver	030 OneGuide KS Fixture Driver	031 OneGuide Path Drill	031 Anchor Drill	031 Mount Driver (OneGuide Anchor Driver)	064 One485 KIT	065 One485 Twist Drill	065 One485 Pilot Drill	066 One485 Drill	070 Denture 4U KIT
032 Anchor Screw	032 OneGuide Twist Drill	032 OneGuide Bone Anchor	033 OneGuide Fixture Anchor	033 CT Checker	071 Denture 4U Guide	071 Posterior Guide	071 Crest Remover	072 Anchor Screw	072 Anchor Drill
033 OneGuide Reamer Drill	044 OneGuide Accessory KIT	046 One Positioning KIT	052 OneMS KIT	053 Tissue Punch	072 Guided Initial Drill	073 Twist Drill	073 Countersink	073 Indicator	073 Path Checker
053 Flattening Drill	054 OneMS Drill	054 OneMS Cortical Drill	055 MS Narrow Ridge NoMount Driver	055 NoMount Driver (TS)	074 Simple Mount Driver	074 Multi Abutment Machine Driver	074 Multi Abutment Outer Driver	076 Positioning Guide KIT	077 Positioning Guide Full KIT

KIT Contents 2/4

078 Guide Drill	078 Single Guide	078 Guide Pin (Fixture)	079 Guide Pin	079 Bridge Guide	112 123 Ultra Twist Drill	113 Parallel Pin (123 Drill)	113 Trial Pin (Ultra-wide)	120 New Hanaro KIT	122 Twist Drill (Stopper Drill)
079 Multi Joint Handle	079 Denture Guide	080 L-wrench	080 Distance Setup Pin	081 Smart Guide KIT	122 Twist Drill (Non stopper Drill)	123 Long Shank Pilot Drill	123 Cortical Drill 2 (TSII, SSII SA)	123 Cortical Drill 3 (Taper Fixture TSIII, SSIII, USIII, KSIII)	123 Countersink (USIII, USII SA, USIII SA Wide PS, Wide)
082 Smart Guide	082 Twist Drill	082 X Sleeve	083 Twist Drill (Ø2.2)	083 Guide Pin	124 Straight Fixture Tap (TSII, USII, SSII SA)	124 Parallel Pin	126 Ultra KIT	127 Direct Drill	127 Cortical Drill (Ultra-wide)
084 122 Taper KIT	085 122 Taper Full KIT	086 122 Taper DrillI	086 Cortical Drill (Ultra-wide)	087 Parallel Pin (122 Taper Drill)	127 Trial Pin (Ultra-wide)	138 485 KIT	139 485 Drill	142 123 Guide Drill	142 Lance Drill (Guide Drill)
094 Taper KIT	095 Taper Ultra KIT	096 Taper Drill	096 Taper Cortical Drill (Taper Fixture TSIII, SSIII, USIII)	097 Taper Ultra Drill	142 Sidecut Drill	143 Drill Extension	143 TS NoMount Driver	143 SS NoMount Driver	144 US NoMount Driver
097 Cortical Drill (Ultra-wide)	097 Parallel Pin (Taper Drill)	098 Tapered Fixture Tap (Taper Fixture TSIII, USIII, SSIII SA)	106 123 Straight Simple KIT	107 123 Twist Drill	144 KS NoMount Driver	144 TS NoMount Torque Driver	145 SS NoMount Torque Driver	145 TS Fixture Driver	145 SS Fixture Driver
107 123 Drill Stopper	108 123 Cortical Drill	110 123 Straight KIT	111 123 Straight Full KIT	112 123 Twist Drill (Stopper Drill)	146 US Fixture Driver	146 KS Fixture Driver	146 Torque Extension	147 Simple Mount Driver	147 Simple Mount Extension

KIT Contents 3/4

147 Simple Open Wrench 	148 Removal Tool (Fixture Mount) 	148 Depth Gauge 	148 Positioning Guide 	148 Tissue Height Gauge (TS) 	161 Osstem Torque Driver 	161 Path Probe (TS) 	161 Path Probe (KS) 	162 Torque Connector 	162 Machine Driver Connector 
149 Ratchet Wrench 	149 L-Wrench 	149 Torque Wrench (Spring Type) 	149 Torque Wrench (Bar Type) 	150 Torque Wrench Set 	162 Driver Handle 	162 Finishing Reamer Set 	163 Reamer Bite 	163 Reamer Tip (Rigid Abutment) 	163 Reamer Tip (Solid, Excellent Solid Abutment) 
150 Tissue Punch 	151 Bone Profiler (TS) 	151 Bone Profiler (US) 	152 Trephine Drill 	152 Machine Driver Handle 	164 CAS KIT 	165 CAS Drill 	165 Guide Drill 	165 Twist Drill (Ø 2.2) 	166 Hydraulic Membrane Lifter Set 
152 Bone Mill 	153 Anterior Hand Driver (Implant) 	153 Torque Handle 	154 Prosthetic Simple KIT 	155 Prosthetic KIT 	166 Stopper 	166 Bone Carrier 	166 Bone Carrier Head 	167 Bone Condenser 	167 Hydraulic Membrane Lifter Tube 
156 Hand Driver 	156 Machine Screw Driver 	157 Torque Driver 	157 Angled Torque Driver 	157 Repair Torque Driver 	167 Membrane Lifter 	168 Depth Gauge 	168 Bone Spreader 	168 Y- Connector 	170 LAS KIT 
158 Solid Abutment Driver 	158 O-ring Abutment Driver 	158 Rigid Outer Driver 	159 Excellent Solid Abutment Driver 	159 Octa Abutment Driver 	171 LAS Full KIT 	172 Dome Drill 	172 Core Drill 	172 Side Wall Drill 	173 Bone Separator 
160 Multi Abutment Machine Driver 	160 Abutment Holder 	160 Abutment Positioning Driver 	160 Multi Abutment Outer Driver 	161 Locator Torque Driver 	173 Stopper 	174 ESSET KIT 	175 Crest Remover 	175 Twist Drill 	175 Saw 



022	OneGuide KIT	1 55	Prosthetic KIT
044	OneGuide Accessory KIT	1 64	CAS KIT
046	OnePositioning KIT	1 70	LAS KIT
052	OneMS KIT	1 71	LAS Full KIT
060	OneCAS KIT	1 74	ESSET KIT
064	One485 KIT	1 78	IM-Cure KIT
070	Denture 4U KIT	1 82	ESR KIT
076	Positioning Guide KIT	1 83	ESR Full KIT
077	Positioning Guide Full KIT	1 92	EFR KIT
081	SmartGuide KIT	1 93	EFR Full KIT
084	122 Taper KIT	1 98	Dr.Cho's Instrument KIT
085	122 Taper Full KIT	1 99	Osstem Basic Instrument KIT
094	Taper KIT	202	Custom KIT
095	Taper Ultra KIT	203	Healing Case
1 06	123 Straight Simple KIT	204	Osteo KIT
1 1 0	123 Straight KIT	205	Osteotome KIT
1 1 1	123 Straight Full KIT	206	Sinus KIT
1 20	New Hanaro KIT	207	Bone Spreader KIT
1 26	Ultra KIT	208	Ridge Split KIT Straight
1 38	485 KIT	209	Ridge Split KIT Offset
1 54	Prosthetic Simple KIT		




Top panel components

Torque Wrench TW30B

Depth Gauge OSDG

For TSIII / IV SSIII USIII KSIII



OneGuide


- Sleeveless type : 2 types, open type and close type
 - Open type can be used in posterior region with limited opening
- Metal sleeve type : 1 close type
 - Inserted to the OneGuide hole for use
 - Option available upon ordering the surgical guide
- 2 guide hole types according to the diameter of the fixture
 - Regular hole (Ø5.1) : F3.5 / 4.0 / 4.5
 - Wide hole (Ø5.8) : F5.0
- Double contact function for excellent implant placement accuracy
 - Drill for double contact with drilling hole and OneGuide
- Simple drilling sequence by using 122 Taper KIT Drill
- Packing unit : surgical guide
 - Option : OneFit Abutment, temporary crown



Tissue Punch RENEWAL 2020

- Used to remove gingiva in flapless surgery
- 7 types according to OneGuide hole diameter
- Other types except the 2 types included in the KIT (OGTP35R, OGTP40W) are sold separately


D	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
Ø3.0	OGTP30R	-
Ø3.5	OGTP35R	-
Ø4.0	OGTP40R	OGTP40WC
Ø4.5	OGTP45R	OGTP45WC
Ø5.0	-	OGTP50WC



Flattening Drill

- Used for narrow or uneven ridges
- Many cutting blades enabling stable removal without bouncing
- 2 types (for below F4.5 / for F5.0)

	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
For below F4.5	OGFD45	-
F5.0	-	OGFD52WC




OneGuide KIT Surgical Instruments

Initial Drill

- Positioning of placement location after using Tissue Punch
- Securing the guide depth of the following drill
- 3 types (F3.5 soft, F4.0/4.5, F5.0)
- Sold separately


	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
F3.5 Soft	OGID2206II	-
F4.0/F4.5	OGID2206I	-
F5.0	-	OGID2206IWC



Initial Drill (Short Type) NEW 2020

- Short type drill with a handle 5.3mm shorter than the Initial Drill
- Used for limited intermaxillary space
- 3 types (F3.5 soft, F4.0/4.5, F5.0)
- Sold separately


	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
F3.5 Soft	OGD2206IIS	-
F4.0/F4.5	OGD2206IS	-
F5.0	-	OGD2206ISWC



OneGuide Twist Drill (Ø2.2)


- Used for placing a F3.5 Fixture in soft bone
- 5 types according to the length

L	TL	Ø2.2
7	36.1	OGTD2207
8.5	36.1	OGTD2208
10	36.1	OGTD2210
11.5	37.6	OGTD2211
13	39.1	OGTD2213



OneGuide Drill

- Taper Drill optimized for III/IV type Fixture
- Used for placing F3.5~F5.0 and 6~13mm Fixture
- Stable drilling with multi-stage structure
- 6mm diameter and F5.5(W) types are sold separately



Regular Hole (Ø5.1)

L	TL	F3.5	F4.0	F4.5
	Y-Dim	0.7	0.9	1.0
6	36.1	OGTPD3506	OGTPD4006	OGTPD4506
7	36.1	OGTPD3507	OGTPD4007	OGTPD4507
8.5	36.1	OGTPD3508	OGTPD4008	OGTPD4508
10	36.1	OGTPD3510	OGTPD4010	OGTPD4510
11.5	37.6	OGTPD3511	OGTPD4011	OGTPD4511
13	39.1	OGTPD3513	OGTPD4013	OGTPD4513

Wide Hole (Ø5.8)


L	TL	F3.5 (W)	F4.5 (W)	F5.0 (W)	F5.5 (W)
	Y-Dim	0.7	1.0	1.0	1.0
6	36.1	OGTPD3506WC	OGTPD4506WC	OGTPD5006WC	OGTPD5506WC
7	36.1	OGTPD3507WC	OGTPD4507WC	OGTPD5007WC	OGTPD5507WC
8.5	36.1	OGTPD3508WC	OGTPD4508WC	OGTPD5008WC	OGTPD5508WC
10	36.1	OGTPD3510WC	OGTPD4510WC	OGTPD5010WC	OGTPD5510WC
11.5	37.6	OGTPD3511WC	OGTPD4511WC	OGTPD5011WC	OGTPD5511WC
13	39.1	OGTPD3513WC	OGTPD4513WC	OGTPD5013WC	OGTPD5513WC

OneGuide KIT Surgical Instruments

RENEWAL 2020

OneGuide Taper Cortical Drill

- Used for placing F4.5 and F5.0 Fixtures in hard bone
- Optimized placement torque by cutting cortical bone
- Product for 13mm diameter Fixtures is sold separately
- Drilling up to the marking line when placing F5.0 6mm



L \ C	Regular Hole (ø5.1)		Wide Hole (ø5.8)	
	Mini	Regular	Mini	Regular
6 / 7 / 8.5mm		OGTCD4507		OGTCD5007WC
10 / 11.5mm		OGTCD4510		OGTCD5010WC
13mm		OGTCD4513		OGTCD5013WC

Fixture Driver


- Used by assembling to a wrench for adjusting the final placement depth
- Yellow groove formed to align the abutment hex direction
- Checked by matching the groove of OneGuide with the groove of driver
- C = Connection



C \	Regular Hole (ø5.1)		Wide Hole (ø5.8)	
	Mini	Regular	Mini	Regular
F3.5	OGFDM50	-	-	-
F4.0 / F4.5	-	OGFDR50	-	-
F5.0	-	-	-	OGFDR57

Fixture Driver (Stopper Type) NEW 2020


- Featuring stopper design to prevent entry below the upper surface of OneGuide hole
- Sold separately
- C = Connection



C \	Regular Hole (ø5.1)		Wide Hole (ø5.8)	
	Mini	Regular	Mini	Regular
F3.5	OGFDSM50	-	-	-
F4.0 / F4.5	-	OGFDSR50	-	-
F5.0	-	-	-	OGFDSR57

NoMount Driver


- Used for placing a NoMount fixture
- ※ It is recommended to place up to 80% of the planned fixture placement depth
- C = Connection



C \	Regular Hole (ø5.1)		Wide Hole (ø5.8)	
	Mini	Regular	Mini	Regular
F3.5	OGNMDM50	-	-	-
F4.0 / F4.5	-	OGNMDR50	-	-
F5.0	-	-	-	OGNMDR57

OneGuide SS NoMount Driver


- Used for placing a SSIII NoMount fixture
- It is recommended to place up to 80% of the planned fixture placement depth
- Sold separately
- P = Platform



P	Regular Hole (ø5.1)	Wide Hole (ø5.8)
	Regular	Regular
For below F4.5	OGNMDR50S	-
F5.0	-	OGNMDR57S

OneGuide US NoMount Driver


- Used for placing a USIII NoMount Fixture
- It is recommended to place up to 80% of the planned fixture placement depth
- Sold separately
- P = Platform



P	Mini	Regular Hole (ø5.1)	Wide Hole (ø5.8)
		Regular	Wide
F3.5	OGNMDM50U	-	-
F4.0 / F4.5	-	OGNMDR50U	-
F5.0	-	-	OGNMDW57U

OneGuide SS Fixture Driver


- Used by assembling to a wrench for adjusting the final placement depth
- Placing SSIII G/H 2.8 up to the octa custom groove marking line
- Yellow groove formed to align the abutment octa direction
- Checked by matching the groove of OneGuide with the groove of driver
- Sold separately
- P = Platform



P	Regular Hole (ø5.1)	Wide Hole (ø5.8)
	Regular	Regular
For below F4.5	OGFDR50S	-
F5.0	-	OGFDR57S

OneGuide US Fixture Driver

- Used by assembling to a wrench for adjustment of the final placement depth
- Yellow groove formed to align the abutment hex direction
- Checked by matching the groove of OneGuide with the groove of driver
- Sold separately
- P = Platform




P	Mini	Regular Hole (ø5.1)	Wide Hole (ø5.8)
		Regular	Wide
F3.5	OGFDM50U	-	-
F4.0 / F4.5	-	OGFDR50U	-
F5.0	-	-	OGFDW57U

OneGuide KIT Surgical Instruments

OneGuide KS NoMount Driver NEW 2020


- Used for placing a KS NoMount Fixture
- It is recommended to place up to 80% of the planned fixture placement depth
- Sold separately
- C = Connection



	Regular Hole (Ø5.1) Regular	Wide Hole (Ø5.8) Regular
C		
For below F4.5	OGNMDR50K	-
F5.0	-	OGNMDR57K

OneGuide KS Fixture Driver


- Used by assembling to a wrench for adjustment of the final placement depth
- Yellow groove formed to align the abutment hex direction
- Checked by matching the groove of OneGuide with the groove of driver
- Below F4.5 : Up to the marking line
 - F3.5 : Up to the lower line, placing the a fixture up to the lower part of the hex custom groove line
- Sold separately
- C = Connection



	Regular Hole (Ø5.1) Regular	Wide Hole (Ø5.8) Regular
C		
For below F4.5	OGFDR50K	-
F5.0	-	OGFDR57K

OneGuide Path Drill 12.2018

- Drill to correct the path deviation during OneGuide surgery
- Drill to form fixture placement path for extraction case
- Flat blade design optimized for cutting inclined bones
- 4 types for each OneGuide hole diameter, 8 types in total : Regular hole (Ø5.1) / Wide hole (Ø5.8)
- Default KIT components : Regular hole (Ø5.1) - Ø2.5 / Wide hole (Ø5.8) - Ø3.0
- 13mm type product adjusts depth according to the marking line (Top line 11.5mm, Midline 10mm, Bottom line 8.5mm)



Regular Hole (Ø5.1)

L \ D	Ø2.5	Ø3.0
7	OGSD2507	OGSD3007
13.0	OGSD2513	OGSD3013

Wide Hole (Ø5.8)

L \ D	Ø2.5	Ø3.0
7	OGSD2507WC	OGSD3007WC
13.0	OGSD2513WC	OGSD3013WC

Anchor Drill

- Used for drilling before using an Anchor Screw

	QGATD13
--	---------



Mount Driver (OneGuide Anchor Driver)

- Used by connecting to a simple mount for placing a fixture (Short type)
- Used by connecting to an Anchor Screw for OneGuide surgery


	ASMDS
--	-------



OneGuide KIT Surgical Instruments

Anchor Screw


- Used for fixing OneGuide in place (e.g. edentulous case)
- Applied selectively in preoperative planning stage



OneGuide Fixture Anchor

11.2019


- Used for fixing OneGuide in place vertically (e.g. edentulous case)
- Placed to the fixture vertically to fix OneGuide in place
- Tightened with 1.2 hex driver (hand mode)
- Only used for a Regular connection of F4.0 or greater
- Sold separately



OneGuide Twist Drill

11.2019

- Used for drilling before using an OneGuide Bone Anchor
- Sold separately



CT Checker

08.2019


- Checking the drilling path through CT scan by connecting to the guide hole before OneGuide procedure (e.g. edentulous case)
- 1 type each for each hole diameter
- Sold separately
- 1 set = 5ea



OneGuide Bone Anchor

11.2019



- Used for fixing OneGuide in place vertically (e.g. edentulous case)
- Mounted on alveolar bone to fix OneGuide in place
- Soft bone : placed directly
- Normal/hard bone : placed after using the OneGuide Twist Drill for Bone Anchor
- Tightened 20rpm FWD with Anchor Driver
- Sold separately



OneGuide Reamer Drill

2019

- Reamer for hole size adjustment after OneGuide template output
- 3 types according to the OneGuide hole size
- Sold separately




Drilling Sequence

OneGuide Drill

TSIII | SSIII | USIII | KSIII


(Length : 10mm)

Ø 3.5




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (Ø2.2)	Drill (F3.5)	Drill (F4.0)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	(F3.5 Soft) ☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒		☒			
Hard	☒	(☒)	☒		☒	☒		

Ø 5.0




Bone Quality	Tissue Punch	Flattening Drill (W)	Initial Drill (W)	Drill (W) (F3.5)	Drill (W) (F4.5)	Drill (W) (F5.0)	Cortical (W) (F5.0)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒		☒			
Hard	☒	(☒)	☒	☒		☒	☒		

Ø 4.0



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒	☒			
Hard	☒	(☒)	☒	☒		☒		

Ø 4.5




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Cortical (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒		☒			
Hard	☒	(☒)	☒	☒		☒	☒		

Drilling Sequence OneGuide Drill


TSIII | SSIII | USIII | KSIII

(Length : 10mm)


G/H 1.8
Ø 3.5




G/H 2.8
Ø 3.5




G/H 1.8
Ø 4.0




G/H 2.8
Ø 4.0



G/H 1.8
Ø 4.5



G/H 2.8
Ø 4.5




Drilling Sequence

OneGuide Drill

TSIII | SSIII | USIII | KSIII


(Length : 10mm)

Ø 3.5




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (Ø2.2)	Drill (F3.5)	Drill (F4.0)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	(F3.5 Soft) ☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒		☒			
Hard	☒	(☒)	☒		☒	☒		

Ø 5.0




Bone Quality	Tissue Punch	Flattening Drill (W)	Initial Drill (W)	Drill (W) (F3.5)	Drill (W) (F4.5)	Drill (W) (F5.0)	Cortical (W) (F5.0)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒		☒			
Hard	☒	(☒)	☒	☒		☒	☒		

Ø 4.0



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒			Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒	☒			
Hard	☒	(☒)	☒	☒		☒		

Ø 4.5



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Cortical (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒				Implant Placement (Up to 80%)	Implant Placement
Normal	☒	(☒)	☒	☒		☒			
Hard	☒	(☒)	☒	☒		☒	☒		


Drilling Sequence

OneGuide Drill


TSIII | SSIII | USIII | KSIII

(Length : 10mm)


Ø 3.5




Ø 5.0



Ø 4.0




Ø 4.5



Drilling Sequence OneGuide Drill


TSIV
(Length : 10mm)

Ø 4.0




Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒				
Normal	☒	(☒)	☒	☒	☒		Implant Placement (Up to 80%)	Implant Placement
Hard								

Ø 4.5



Bone Quality	Tissue Punch	Flattening Drill	Initial Drill	Drill (F3.5)	Drill (F4.0)	Drill (F4.5)	Cortical (F4.5)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒	☒				
Normal	☒	(☒)	☒	☒		☒		Implant Placement (Up to 80%)	Implant Placement
Hard									

Ø 5.0



Bone Quality	Tissue Punch	Flattening Drill (W)	Initial Drill (W)	Drill (W) (F3.5)	Drill (W) (F4.5)	Drill (W) (F5.0)	Cortical (W) (F5.0)	Nomount Driver	Fixture Driver
Soft	☒	(☒)	☒	☒	☒				
Normal	☒	(☒)	☒	☒		☒		Implant Placement (Up to 80%)	Implant Placement
Hard									


OSSTEM[®]
IMPLANT

OneGuide Accessory KIT

(OOGAK)

NEW 2020

- KIT consisted of the tools selected by user
- Possible to accommodate the products not included in the OneGuide KIT by default such as OneGuide Bone/Fixture Anchor, and SS/US/KS Driver
- Spare holes deployed by rubber size (Large 4, Medium 8, Small 4) for user preferences




OSSTEM[®]
IMPLANT

OnePositioning KIT (OOPK) 2019

For TSIII / IV SSIII USIII / IV KSIII III / IV Ultra-wide

- Selecting the initial placement position, path and depth using OneGuide
- Removing OneGuide after F3.5 drilling and proceeding up to fixture placement through manual surgery



OSSTEM[®]
IMPLANT

Drilling Sequence

OneGuide Drill

+ 122 Taper Drill

TSIII/IV

SSIII

USIII/IV

KSIII

III/IV Ultra-wide

(Length : 10mm)

Ø 3.5

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch	Flattening	Initial	F3.5	Guide Drill	F4.0	
Soft	(☒)	(☒)	☒ F3.5 soft (Option)		☒		Implant Placement
Normal	(☒)	(☒)	☒	☒			
Hard	(☒)	(☒)	☒	☒		☒	

Ø 4.0

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch	Flattening	Initial	F3.5	F4.0	F4.5	
Soft	(☒)	(☒)	☒	☒			Implant Placement
Normal	(☒)	(☒)	☒	☒	☒		
Hard	(☒)	(☒)	☒	☒		☒	

Ø 4.5

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch	Flattening	Initial	F3.5	F4.0	F4.5	F4.5	
Soft	(☒)	(☒)	☒	☒	☒			Implant Placement
Normal	(☒)	(☒)	☒	☒		☒		
Hard	(☒)	(☒)	☒	☒			☒	

Ø 5.0

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch (W)	Flattening (W)	Initial (W)	F3.5 (W)	F4.5	F5.0	F5.5	
Soft	(☒)	(☒)	☒	☒	☒			Implant Placement
Normal	(☒)	(☒)	☒	☒		☒		
Hard	(☒)	(☒)	☒	☒			☒	

Ø 5.5

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch (W)	Flattening (W)	Initial (W)	F3.5 (W)	F5.0	F5.5	F5.5 Cortical	
Soft	(☒)	(☒)	☒	☒	☒			Implant Placement
Normal	(☒)	(☒)	☒	☒		☒		
Hard	(☒)	(☒)	☒	☒		☒	☒	

Ø 6.0

OnePositioning KIT

122 Taper KIT

Bone Quality	Tissue Punch (W)	Flattening (W)	Initial (W)	F3.5 (W)	F5.0	F5.5	F6.0	F6.0 Cortical	
Soft	(☒)	(☒)	☒	☒	☒	☒			Implant Placement
Normal	(☒)	(☒)	☒	☒	☒		☒		
Hard	(☒)	(☒)	☒	☒	☒		☒	☒	

Drilling Sequence

OneGuide Drill

+ 122 Taper Drill

TSIII/IV

SSIII

USIII/IV

KSIII

III/IV Ultra-wide

(Length : 10mm)

Ø 7.0

OnePositioning KIT

122 Taper KIT

*Up to Ø5.0 with OneGuide planning

Bone Density	Tissue Punch (W)	Flattening (W)	Initial (W)	F3.5 (W)	F5.0	F6.0	F7.0	F7.0 Cortical	Implant Placement
Soft	(☒)	(☒)	☒	☒	☒	☒			
Normal	(☒)	(☒)	☒	☒	☒		☒		
Hard	(☒)	(☒)	☒	☒	☒		☒	☒	



OSSTEM[®]
IMPLANT

For **TSIII** **MS**
Ø3.0

Top panel components

Torque Wrench
TW30B

Depth Gauge
MSDG





OneGuide Template

- Sleeveless type : 2 types, open type and close type
 - Open type can be used in posterior region with limited opening
- Metal sleeve type : 1 close type
 - Placed to the OneGuide hole for use
 - Option available upon ordering the surgical guide
- 1 guide hole type for narrow fixture diameter
 - Narrow hole (Ø3.6) : MS narrow Ø2.0 / 2.5 / 3.0, TSIII Ø3.0
- Double contact function for excellent positioning accuracy
 - Drill for double contact with drilling hole and OneGuide
- Simple drilling sequence by using conventional drilling sequence
- Packing unit : surgical guide
 - Option : temporary crown





Tissue Punch **RENEWAL** 2020

- Used to remove gingiva in flapless surgery




Flattening Drill

- Used for narrow or uneven ridges
- Many cutting edges enabling stable removal without bouncing



OneMS Drill

- Optimized Straight Drill for MS implant / TSIII Ø3.0 Fixture
(For placing MS Ø2.0~3.0, TSIII Ø3.0 Fixtures)
- OneMS Cortical Drill used for placing a TSIII Ø3.0 Fixture in hard bone
- Recommend using 8.5mm Drill within the same diameter for stable drilling
(Inducing the double contact feature)





Narrow Hole (Ø3.6)

L	TL	Ø1.5	Ø1.8	Ø2.3	Ø2.7
8.5	37.5	OMD1508	OMD1808	OMD2308	OMD2708
10	39.0	OMD1510	OMD1810	OMD2310	OMD2710
11.5	40.5	OMD1511	OMD1811	OMD2311	OMD2711
13	42.0	OMD1513	OMD1813	OMD2313	OMD2713


OneMS Cortical Drill

- Drill used for removing cortical bone from hard bone
- Drill used for expanding the cortical bone after using the Straight Drill (for TSIII Ø3.0 Fixture only)



MS Narrow Ridge NoMount Driver

- Used for placing a MS implant Narrow Ridge
- Used by matching the triangular marking with the side of the implant



NoMount Driver (TS)

- Used for placing a TSIII Ø3.0 NoMount Fixture
- It is recommended to insert up to 80% of the planned fixture placement depth
- C = Connection



MS Narrow Ridge Fixture Driver

- Used by assembling to a wrench for adjusting the final placement depth of a MS implant Narrow Ridge
- Used by matching the triangular marking with the side of the implant
- Placing up to the lower marking line for G/H 4.0



OneMS KIT Surgical Instruments

Fixture Driver (TS)


- Used by assembling to a wrench for adjusting the final placement depth of a TSIII Ø3.0 Fixture
- Yellow groove formed to align the abutment hex direction
- Checked by matching the groove of OneGuide with the groove of driver
- C = Connection

C		Narrow Hole (Ø3.6)
		Mini
TSIII Ø3.0		OGFDM35



Driver Separator

- When the driver is caught after MS Implant placement, insert the driver separator into the driver groove and remove it by using the lever principle



Fixture Driver (TS, Stopper Type) NEW 2020

- Featuring stopper design to prevent entry below the upper surface of OneGuide hole
- Sold separately
- C = Connection

C		Narrow Hole (Ø3.6)
		Mini
TSIII Ø3.0		OGFDSM35




OneMS Path Drill 12.2018

- Drill to correct the path deviation during OneGuide surgery
- Drill to form fixture placement path for extraction case
- Flat blade design optimized for cutting inclined bones
- 2 types for each drill diameter, 4 types in total : Narrow hole (Ø3.6)
- 13mm type product adjusts depth according to the marking line (Top line 11.5mm, Midline 10mm, Bottom line 8.5mm)

Narrow Hole (Ø3.6)

L \ D	Ø1.5	Ø2.0
7.0	OMSD1507	OMSD2007
13.0	OMSD1513	OMSD2013



Adapter

- Driver enabling the driver for Torque Wrench to be used for engine




OneMS Lance Drill

- Forming a hole in bone to facilitate initial drilling
- Bone density can be checked through drilling
- Sold separately

Narrow Hole (Ø3.6)

L \ TL	Ø1.5
8.5	37.5 OMLD1508
10	39.0 OMLD1510
11.5	40.5 OMLD1511
13	42.0 OMLD1513



Drilling Sequence OneMS Drill


TSIII (Ø3.0) | MS
(Length : 10mm)

MS Ø 2.0



Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.5)	NoMount Driver	Fixture Driver		
					G/H 2.5	G/H 4.0	Denture
Soft	☑	(☒)	☑	☑		☑	
Normal	☑	(☒)	☑	☑		☑	
Hard	☑	(☒)	☑	☑		☑	

TSIII Ø 3.0




Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø2.3)	OneMS Drill (Ø2.7)	F4.5 Cortical Drill	NoMount Driver	Fixture Driver
Soft	☑	(☒)	☑			☑	☑
Normal	☑	(☒)	☑	☑		☑	☑
Hard	☑	(☒)	☑	☑	☑	☑	☑

MS Ø 2.5



Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	NoMount Driver	Fixture Driver		
						G/H 2.5	G/H 4.0	Denture
Soft	☑	(☒)	☑	-	☑		☑	
Normal	☑	(☒)	☑	-	☑		☑	
Hard	☑	(☒)	-	☑	☑		☑	

MS Ø 3.0



Bone Quality	Tissue Punch	Flattening Drill	OneMS Drill (Ø1.8)	OneMS Drill (Ø2.3)	OneMS Drill (Ø2.7)	NoMount Driver	Fixture Driver		
							G/H 2.5	G/H 4.0	Denture
Soft	☑		☑	-	-	☑		☑	
Normal	☑		☑	-	-	☑		☑	
Hard	☑		-	☑	☑	☑		☑	

In case of Fixture 10 / 11.5 / 13mm sequence, precede number of drilling with 8.5mm Drill for each step

Ex. Ø2.5×11.5mm MS Fixture

: Tissue Punch ► Flattening Drill ► Ø1.8×8.5mm ► Ø1.8×11.5mm ► NoMount Driver ► Fixture Driver

For

TSII / III

SSII / III

USII / III

KSIII

Lower panel components

Bone Carrier Head
OCBCH32, OCBCH37W

Depth Gauge
OCDG

Depth Gauge (W)
OCDGW

Bone Carrier
OCBC530

Hydraulic Membrane Lifter Tube
SNMT

Bone Condenser
SNBC1114

Top panel components

OGFD45

OCD2207
OCD2210

OCD2807
OCD2810

OCD3107
OCD3110

OCD3307
OCD3310

OCD3607
OCD3610

Flattening Drill

Twist Drill (Ø2.2)

OneCAS Drill (Ø2.8)

OneCAS Drill (Ø3.1)

OneCAS Drill (Ø3.3)

OneCAS Drill (Ø3.6)

Stopper
OCDS01
OCDS04
OCDS07
Yellow

Stopper
OCDS02
OCDS05
OCDS08
Purple

Stopper
OCDS03
OCDS06
OCDS09
Blue

Stopper(W)
OCDS02W
OCDS05W
OCDS08W
Purple

Stopper(W)
OCDS01W
OCDS04W
OCDS07W
Yellow

Stopper(W)
OCDS03W
OCDS06W
OCDS09W
Blue

OneCAS Hydraulic Membrane Lifter
OCHML

Flattening Drill (W)
OGFD52WC

Twist Drill (W) (Ø2.2)
OCD2207WC
OCD2210WC

OneCAS Drill (W) (Ø2.8)
OCD2807WC
OCD2810WC

OneCAS Drill (W) (Ø3.1)
OCD3107WC
OCD3110WC


OneCAS Hydraulic Membrane Lifter(W)
OCHMLW

OneCAS Drill (W) (Ø3.8)
OCD3807WC
OCD3810WC

OneCAS Drill (W) (Ø4.1)
OCD4107WC
OCD4110WC

OneCAS Twist Drill (Ø2.2)

- Drilling 1mm under the depth to maxillary sinus floor
- Used with a stopper for safe lifting
- 1mm shorter than a normal Twist Drill



Regular Hole (Ø5.1)


L	TL	Ø2.2
	Y-Dim	0.6
7	33.2	OCD2207
10	36.2	OCD2210

Wide Hole (Ø5.8)

L	TL	Ø2.2
	Y-Dim	0.6
7	33.2	OCD2207WC
10	36.2	OCD2210WC

OneCAS Drill

- Used with a guide of OneGuide system
- Safe lifting of the membrane for maxillary sinus procedure
- Used at low speed for autogenous bone collection
- Used with a stopper for safe lifting
- Final drill diameter selected based on the bone quality
- Recommended speed : 400~800rpm
- 4 type drills of Ø3.3 and Ø3.6 sold separately (OCD3307WC, OCD3310WC, OCD3607WC, OCD3610WC)



Regular Hole (Ø5.1)

L	TL	Ø2.8	Ø3.1	Ø3.3	Ø3.6
7	33.6	OCD2807	OCD3107	OCD3307	OCD3607
10	36.6	OCD2810	OCD3110	OCD3310	OCD3610

Wide Hole (Ø5.8)

L	TL	Ø2.8 (W)	Ø3.1 (W)	Ø3.3 (W)	Ø3.6 (W)	Ø3.8 (W)	Ø4.1 (W)
7	33.6	OCD2807WC	OCD3107WC	OCD3307WC	OCD3607WC	OCD3807WC	OCD4107WC
10	36.6	OCD2810WC	OCD3110WC	OCD3310WC	OCD3610WC	OCD3810WC	OCD4110WC

OneCAS KIT Surgical Instruments

OneCAS Stopper

- Number marking on the stopper indicates the stopping distance for drilling or tool assembly
- Check in the mid panel of the kit, protruding length marked in blue at connecting 7mm drill and protruding length marked in red at connecting 10mm drill
- Apply color coding by length
- Recommended use cycle : 50 times

Regular Hole (Ø5.1)

	L	1	2	3	4	5	6	7	8	9
										
		OCDS01	OCDS02	OCDS03	OCDS04	OCDS05	OCDS06	OCDS07	OCDS08	OCDS09
Color		Yellow	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue

Wide Hole (Ø5.8)

	L	1	2	3	4	5	6	7	8	9
										
		OCDS01W	OCDS02W	OCDS03W	OCDS04W	OCDS05W	OCDS06W	OCDS07W	OCDS08W	OCDS09W
Color		Yellow	Purple	Blue	Yellow	Purple	Blue	Yellow	Purple	Blue

Depth Gauge

- Checking the internal lifting of sinus
- Measuring residual bone depth
- Used with a stopper for safe lifting
- Marking line of the same depth as 10mm drill




	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
	OCDG	OCDGW

Hydraulic Membrane Lifter

NEW 2020

- Dedicated maxillary sinus hydraulic lifting instrument for OneCAS KIT
- Hydraulic pressure is used to separate and lift the sinus membrane
- Used by placing the body until the marking line meets the upper surface of OneGuide hole
- Winged design with optimized sealing for flapless procedure

	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
	OCHML	OCHMLW




Bone Carrier Head

NEW 2020

- Dedicated maxillary sinus filling instrument for OneCAS KIT
- Used by placing into the OneGuide hole to the end
- OCBCH30 : Used after drilling with OneCAS Drill Ø3.1
- OCBCH32 : Used after drilling with OneCAS Drill Ø3.3/Ø3.6
- OCBCH37W : Used after drilling with OneCAS Drill Ø3.8/Ø4.1
- Used repeatedly by filling bone material in the back of the marking line of the head and taking little by little with a bone condenser to completely fill the inside of the maxillary sinus

	Regular Hole (Ø5.1)	Wide Hole (Ø5.8)
	OCBCH30 OCBCH32	OCBCH37W



Bone Carrier

NEW 2020

- Dedicated maxillary sinus filling instrument for OneCAS KIT
- Mounting the head by fastening the handle in the back of the body
- Replaceable head for use

	OCBCS30
--	---------




For TSIII SSIII USIII KSIII

Top panel components


Torque Wrench
TW30B

Depth Gauge
OSDG



One485 Twist Drill

- Initial Drill for determining the placement position and ensuring the guide depth for other drills
- Drilling with a straight blade up to -1mm of the fixture placement depth
- 5 types according to the OneGuide hole diameter
 - Regular hole (Ø5.1) / Wide hole (Ø5.8)
- F4.0/4.5 6mm (extra short implant) type are used 7mm drill
 - Bottom marking line : 6mm, Top line : 7mm
- F5.0 4mm, 5mm (extra short implant) type are used 6mm drill
 - Bottom marking line : 4mm, Midline : 5mm, Top line : 6mm
- Recommended speed : 800 or 1,200rpm



L	TL	Regular Hole (Ø5.1) F4.0 / F4.5	Wide Hole (Ø5.8) F5.0
6.0	32.4	-	O4D2206W
7.0	33.2	O4D2207	O4D2207W
8.5	34.7	O4D2208	O4D2208W

One485 Pilot Drill

- Medium drill for expanding hole diameter
- Tip blade in the shape of 485 Drill, and the side blade in the shape of tapered drill
- 3 types according to the OneGuide hole diameter
 - Regular hole (Ø5.1) / Wide hole (Ø5.8)
- 4mm drill used for 4-5mm Fixtures, and 6mm drill used for 6~8.5mm Fixtures
- Recommended speed : 800 or 1,200rpm



L	TL	Regular Hole (Ø5.1) F4.5	TL	Wide Hole (Ø5.8) F5.0W
4.0	-	-	33.1	O4D3504W
6.0	33.9	O4D3506	32.9	O4D3506W

One485 KIT Surgical Instruments

One485 Drill

- Final drill for final expansion and placement torque optimization
- Tip blade in the shape of 485 Drill, and the side blade in the shape of tapered drill
- 19 types according to the OneGuide hole diameter
 - Regular hole (Ø5.1) / Wide hole (Ø5.8)
- F4.5 and F5.0 hard drill used for placing F4.5 and F5.0 Fixtures in hard bone
- Recommended speed : 800 or 1,200rpm



Regular Hole (Ø5.1)

L \ TL		F4.0	F4.5	F4.5 Hard
6.0	33.9	O4D4006	O4D4506	O4HD4506
7.0	33.9	O4D4007	O4D4507	O4HD4507
8.5	35.4	O4D4008	O4D4508	O4HD4508

Wide Hole (Ø5.8)

L \ TL		F5.0 (W)	F5.0 (W) Hard
4.0	33.1	O4D5004W	O4HD5004W
5.0	33.1	O4D5005W	O4HD5005W
6.0	32.9	O4D5006W	O4HD5006W
7.0	33.9	O4D5007W	O4HD5007W
8.5	35.4	O4D5008W	O4HD5008W



Drilling Sequence One485 Drill

TSIII


SSIII

USIII

KSIII


(Length : 7mm)

Ø 4.0




Bone Quality	Tissue Punch	Flattening Drill	Twist Drill	Pilot Drill	Drill (F4.0)	Drill (F4.5)	Fixture
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Implant Placement
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	

Ø 4.5



Bone Quality	Tissue Punch	Flattening Drill	Twist Drill	Pilot Drill	Drill (F4.5)	Hard Drill (F4.5)	Fixture
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Implant Placement
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	

Ø 5.0




Bone Quality	Tissue Punch (W)	Flattening Drill (W)	Twist Drill (W)	Pilot Drill (W)	Drill (W) (F5.0)	Hard Drill (W) (F5.0)	Fixture
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Implant Placement
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	

OSSTEM[®]
IMPLANT

Denture 4U KIT (OD4UK) NEW 2020

For TSII / III USII / III



Denture 4U KIT Surgical Instruments

Denture 4U Guide


- Guide for stable and accurate initial and intermediate drilling for Denture 4U procedure
- Anterior guide : drilling positioning for Ø2.2 in anterior region (tooth number 2 and 3 positions marked)
- Posterior guide : drilling positioning for Ø3.0 drill in posterior region
 - ※ Used by assembling with the posterior guide of desired angle
- Removable Denture 4U Guide handle



D4UG

Posterior Guide

- Used by assembling the anterior guide prior to procedure
 - ※ Assembled with the angle marking side shown
- Adjusting the fixture placement position in posterior region and buccolingual inclination angle
- Selecting the angle of the posterior guide through CT scan recommended prior to procedure
 - ※ Replaceable during procedure
- Drilling by slowly entering the guide hole, referring to the marking line on the side of the posterior guide hole
- Drilling depth adjusted by drilling to the bottom marking line in the mesial direction
- Marking line spacing on the rod : 2mm




Marking Bottom Line, Check Mesial Direction

Degree	0°	17°	30°
	D4UPG0	D4UPG17	D4UPG30

Crest Remover

- Used for bone flattening for Denture 4U Guide procedure
- Marking the fixture placement position after removing narrowed ridge
- Recommended speed
 - Angled type : 1,200~1,500rpm
 - Straight type : 15,000~30,000rpm

L	D	Ø5.0
29		CERM50A
45		CERM50S



Anchor Screw

- Used to fix the bone in place by connecting it to the fixed center hole of the Denture 4U Guide and the fixed hole of the posterior guide
- Fixing the Anchor Screw with the Mount Driver; if the Anchor Screw is not fixed well at this time, it should be drilled first using an Anchor Drill
 - ※ Anchor drill used first for normal/hard bone
- Selecting an Anchor Screw of appropriate length according to the degree of posterior bone retraction
- Engine stop to prevent Anchor Screw from spinning with no traction when in contact with the guide

L \ D	Ø1.65
5	D4UAS5
8	D4UAS8
11	D4UAS11



Anchor Drill

- Used to form a hole in normal/hard bone prior to tightening an Anchor Screw
- Drilling with 3mm drill prior to additional drilling with 11mm drill recommended

L \ D	Ø1.65
3	D4UAD3
11	D4UAD11



Guided Initial Drill

- Used for drilling in anterior region : Ø2.2 drilling into the anterior guide hole of the Denture 4U Guide
- Drilling by selecting a desired drilling hole of the anterior guide
- Recommended speed : 800rpm

L \ D	Ø2.2
5	GD2208NC



Twist Drill

- Drilling by slowly entering the guide hole, with the angle matched as much as possible, referring to the marking line on the side of the posterior guide hole
- Drilling depth adjusted by drilling to the bottom marking line in the mesial direction
- Marking line spacing on the rod : 2mm
- Recommended speed : 800rpm

L \ D	Ø3.0
5	D4U2D3005
18	D4U2D3018



Countersink

- Drill for using the Taper Drill after removing the Denture 4U Guide
 - ※ For removing bone interference from the stopper of the Taper Drill
- Removing bone interference upon mounting a Multi Angled Abutment


L \ D	Ø3.0
5	D4UCS



Indicator

- Checking the location of the mental foramen, and the placement direction and length of the fixture beforehand for stable procedure
 - ※ For checking the location of the mental foramen by opening a flap completely

L \ D	Ø3.0
5	D4UI



Path Checker

- Checking the location of the mental foramen by predicting the extended line of the path checker through panoramic or CT scan
 - ※ For checking the location of the mental foramen without opening a flap completely

L \ D	Ø3.0
5	D4UPC



Denture 4U KIT

Surgical Instruments

Simple Mount Driver

- Used for placing an Anchor Screw to stably fix the Denture 4U Guide in place



Multi Abutment Machine Driver

- Dedicated Machine Driver for a Multi Abutment



Multi Abutment Outer Driver

- Dedicated Torque Driver for a Multi Abutment



Positioning Guide KIT

(OPGPK) 07.2015






Positioning Guide Full KIT

(OPGAK) 07.2015



Lower panel components



Positioning Guide KIT


Surgical Instruments

Guide Drill

07.2015

- Initial Drill : For initial drilling, assembled to the single guide to adjust the drilling depth
- Ø2.2 Twist Drill : Used with the bridge guide for initial drilling
- Ø3.0 Twist Drill : For subsequent drilling of Ø2.2 Twist Drill, drilling path guide


	Ø2.2	Ø3.0
Initial Drill	GD2208NC	-
Twist Drill	GD2213FNC	2D3011LC01



Single Guide

07.2015

- Transparent material applied to facilitate the viewing of the position and direction for drilling
- 6 types considering mesio-distal crown diameters (Ø5.0~10.0)
- Packing unit : 2ea
- ※ Drilling depth adjusted to 6, 7 or 8mm using the marking line of the Initial Drill, based on the top line of the single guide
- ※ Disposable, Do not reuse



	F5.0	F6.0	F7.0	F8.0	F9.0	F10.0
	SGB050	SGB060	SGB070	SGB080	SGB090	SGB100

Guide Pin (Fixture)

07.2015

- Pin for checking the path and fixing the single guide in place after placing a fixture
- C = Connection

	Mini	Regular
	PGSPM	PGSPR



Guide Pin

07.2015

- Pin for checking the drilling path and fixing the single guide in place



Bridge Guide

07.2015

- Guide for adjusting the direction and distance for drilling
- Fan type : Selectable in 0.5mm increments (7~12.5mm)
- Compass type : Adjustable in 1 mm increments (5~24mm)
- Used after adjusting the distance in the distance setup of the mid panel of KIT



Type \ Distance	7~8.5	9~10.5	11~12.5	5~24
Fan	PGBRA070	PGBRA090	PGBRA110	-
Compass	-	-	-	PGBPA

Multi Joint Handle

Option 07.2015

- Instrument to place the guide from the outside of the oral cavity by connecting to the ball head of the bridge guide



Denture Guide

Option 07.2015

- Guide with adjustable angle for each patient in edentulous case
- Drilling in the oral cavity with the angle fixed with an L-wrench in working model model after adjusting the angle according to the arch shape of the patient
- Marking line refers to the No. 2,3,4,5,6 positions from the center



Positioning Guide KIT

Surgical Instruments

L-wrench Option 07.2015

- Instrument to adjust the size of the denture guide and keep it in place



LWC20H

Distance Setup Pin Option 07.2015

- Pin for bridge guide compass type and denture guide fixation



PGDSP

SmartGuide KIT

(OSGK) 12.2015



Lower panel components

Guide Pin (4ea)
SGP22



Round bur (2ea)
RAHM1018




Twist Drill (2ea) (Ø2.2) or stone models
2D2208LC01



SGCB30S

X Sleeve
(4ea)

Twist Drill
(Ø2.2)
SGTD2207S



Twist Drill
(Ø3.0)
QGTD3008



SmartGuide KIT Surgical Instruments

SmartGuide 12.2015

- Thermoplastic surgical guide
 - Freely deformable after immersion in about 70°C water for about 1 minute
 - Curing at room temperature after 1 minute from deformation
- ※ Disposable, Do not reuse, Use after low temperature disinfection (Do not autoclave or hydrogen peroxide)

Type	Single	Free-end Bridge	2-Unit Br.: small	2-Unit Br.: large
				
	SGTSS	SGTFB90LS	SGTB63SS	SGTB85LS

Twist Drill (Ø2.2) For stone models 12.2015

- Used for initial marking on the working model
- Use cycle : 10 times
- Additional drilling after using the round bur
- Recommended speed : 1,200~1,500rpm

D	Ø2.2
	2D2208LC01



Twist Drill 12.2015

- Drill used through the guide in the oral cavity
- Stable drilling by connecting to the SmartGuide sleeve
- After initial drilling with Ø2.2 drill, additional drilling with Ø3.0 drill
- Recommended speed : 1,200~1,500rpm

D	Ø2.2	Ø3.0
	SGTD2207S	QGTD3008



Guide Pin 12.2015

- Assembled to the working model for fixing the SmartGuide in place
- Connected to the SmartGuide sleeve

	SGP22
--	-------



X Sleeve 12.2015

- Instrument to check if th guide is produced as intended through CT scan or x-ray by connecting to the SmartGuide sleeve
- After connecting to the SmartGuide outside the oral cavity, assemble in the oral cavity

	SGCB30S
--	---------





122 Taper KIT (0122TPK) 09.2016



Top panel components



For TSIII / IV SSIII USIII / IV KSIII



122 Taper Full KIT (0122TPFK) 01.2018

Top panel components

For TSIII / IV SSIII USIII / IV KSIII III / IV Ultra-wide




122 Taper KIT

Surgical Instruments


122 Taper Drill

- Dedicated Taper Drill for Taper (III type) Fixture
- Types available for each diameter and length
- Color coded handle indicating the fixture diameter
- Drill slightly larger in diameter used for removing cortical bone from hard bone
- Included in 122 Taper KIT only (not included in Taper KIT)
- F = Fixture



Parallel Pin (122 Taper Drill)

- Dedicated Parallel Pin for 122 Taper Drill
- Used for checking the position and direction of bone preparation
- Bottom part for Ø2.2 drill, and top part for guide drill
- Included in 122 Taper KIT only (not included in Taper KIT)
- Other components same as Taper KIT



APP2227


※ Refer to surgical instruments for other components (from p142)

L	TL	F3.5	F4.0	F4.5	F5.0	F5.5	F6.0	F7.0
	Y-Dim.	0.7	0.9	1.0	1.0	1.0	1.0	1.0
4.0	29.5	122TPD 3504	122TPD 4004	122TPD 4504	122TPD 5004	122TPD 5504	-	-
5.0	29.5	122TPD 3505	122TPD 4005	122TPD 4505	122TPD 5005	122TPD 5505	-	-
6.0	30.5	122TPD 3506	122TPD 4006	122TPD 4506	122TPD 5006	122TPD 5506	122TPD 6006	122TPD 7006
7.0	31.5	122TPD 3507	122TPD 4007	122TPD 4507	122TPD 5007	122TPD 5507	122TPD 6007	122TPD 7007
8.5	33	122TPD 3508	122TPD 4008	122TPD 4508	122TPD 5008	122TPD 5508	122TPD 6008	122TPD 7008
10	34.5	122TPD 3510	122TPD 4010	122TPD 4510	122TPD 5010	122TPD 5510	122TPD 6010	122TPD 7010
11.5	34.5	122TPD 3511	122TPD 4011	122TPD 4511	122TPD 5011	122TPD 5511	122TPD 6011	122TPD 7011
13	36	122TPD 3513	122TPD 4013	122TPD 4513	122TPD 5013	122TPD 5513	122TPD 6013	122TPD 7013
15	38	122TPD 3515	122TPD 4015	122TPD 4515	122TPD 5015	122TPD 5515	-	-
Color		Yellow	Green	Blue	Red	Yellow	Green	Blue

Cortical Drill (Ultra-wide) 01.2009

- Drill used for removing cortical bone from hard bone (for Ultra-wide)
- Dedicated drill for each fixture diameter
- Drilling up to the lower marking line recommended
- F = Fixture

F6.0	F7.0
CD4C60	CD4C70




Drilling Sequence 122 Taper Drill


TSIII | SSIII | USIII | KSIII

(Length : 10mm)


Ø 3.0mm



Ø 3.5mm



Ø 4.0mm




F5.5 Taper Cortical Drill marking bottom line for 6mm Fixtures, midline for 7mm Fixtures, top line for 8.5mm or greater Fixtures


Recommended placement torque ≤ 40Ncm

TS Fixture placed to a depth 1mm deeper than the bone level for normal bone, and to the bone level for soft bone to maintain fixation stability


Ø 4.5mm



Ø 5.0mm



Ø 5.5mm




Drilling Sequence 122 Taper Drill

TSIII Ultra-wide | SSIII Ultra-wide | USIII Ultra-wide


(Length : 10mm)

Ø 6.0mm



Bone Quality	Guide Drill	122 Taper Drill (F3.5)	122 Taper Drill (F5.0)	122 Taper Drill (F5.5)	122 Taper Drill (F6.0)	Taper Cortical Drill (F6.0)	Ø6.0 Fixture
Soft	☒		☒	☒			Implant Placement
Normal	☒	☒	☒		☒		
Hard	☒	☒	☒		☒	☒	

Ø 7.0mm



Bone Quality	Guide Drill	122 Taper Drill (F3.5)	122 Taper Drill (F5.0)	122 Taper Drill (F6.0)	122 Taper Drill (F7.0)	Taper Cortical Drill (F7.0)	Ø7.0 Fixture
Soft	☒		☒	☒			Implant Placement
Normal	☒	☒	☒		☒		
Hard	☒	☒	☒		☒	☒	

Drilling Sequence 122 Taper Drill

TSIV | USIV


(Length : 10mm)

Ø 4.0mm




Bone Quality	Twist Drill (Ø2.2)	122 Taper Drill (F3.5)	Ø4.0 Fixture
D4	☒		Implant Placement
Soft	☒	☒	

Ø 4.5mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	122 Taper Drill (F4.0)	Ø4.5 Fixture
D4		☒		Implant Placement
Soft	☒		☒	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	122 Taper Drill (F4.5)	Ø5.0 Fixture
D4		☒		Implant Placement
Soft	☒		☒	


F5.5 Taper Cortical Drill marking bottom line for 6mm Fixtures, midline for 7mm Fixtures, top line for 8.5mm or greater Fixtures
Recommended placement torque ≤ 40Ncm
TS Fixture placed to a depth 1mm deeper than the bone level for normal bone, and to the bone level for soft bone to maintain fixation stability

Drilling Sequence 122 Taper Drill

TSIV Ultra-wide


(Length : 10mm)

Ø 6.0mm



Bone Quality	Guide Drill	122 Taper Drill (F5.0)	122 Taper Drill (F5.5)	Ø6.0 Fixture
D4	☒	☒		Implant Placement
Soft	☒	☒	☒	

Ø 7.0mm



Bone Quality	Guide Drill	122 Taper Drill (F5.0)	122 Taper Drill (F5.5)	122 Taper Drill (F6.0)	Ø7.0 Fixture
D4	☒	☒	☒		Implant Placement
Soft	☒	☒	☒	☒	



F5.5 Taper Cortical Drill marking bottom line for 6mm Fixtures, midline for 7mm Fixtures, top line for 8.5mm or greater Fixtures

Recommended placement torque ≤ 40Ncm



TS Fixture placed to a depth 1mm deeper than the bone level for normal bone, and to the bone level for soft bone to maintain fixation stability

Taper KIT (OTSK) 09.2016



Top panel components



For TSIII / IV SSIII USIII / IV KSIII



Taper Ultra KIT (HULTPK) 07.2013

Lower panel components


For III Ultra-wide



Taper KIT Surgical Instruments

Taper Drill

- Dedicated Taper Drill for Taper(III type) fixtures of each diameter and length
- Stopper Drill with 1mm margin
- Color coded handle indicating the fixture diameter
- F3.5 : Yellow, F4.0 : Green, F4.5 : Blue, F5.0 : Red, F5.5 : Yellow
- Included in Taper KIT only (not included in 122 Taper KIT)



L	TL	F3.5	F4.0	F4.5	F5.0	F5.5
	Y-Dim.	0.8	0.9	1.0	1.0	1.0
5.0	29.5	TPD3C 3505	TPD3C 4005	TPD3C 4505	TPD3C 5005	-
6.0	30.5	TPD3C 3506	TPD3C 4006	TPD3C 4506	TPD3C 5006	TPD3C 5506
7.0	31.5	TPD3C 3507	TPD3C 4007	TPD3C 4507	TPD3C 5007	TPD3C 5507
8.5	33	TPD3C 3508	TPD3C 4008	TPD3C 4508	TPD3C 5008	TPD3C 5508
10	34.5	TPD3C 3510	TPD3C 4010	TPD3C 4510	TPD3C 5010	TPD3C 5510
11.5	34.5	TPD3C 3511	TPD3C 4011	TPD3C 4511	TPD3C 5011	TPD3C 5511
13	36	TPD3C 3513	TPD3C 4013	TPD3C 4513	TPD3C 5013	TPD3C 5513
15	38	TPD3C 3515	TPD3C 4015	TPD3C 4515	TPD3C 5015	TPD3C 5515
Color		Yellow	Green	Blue	Red	Yellow

Taper Cortical Drill
(Taper Fixture TSIII, SSIII, USIII)

- Drill used for removing cortical bone from hard bone (used right after Taper Drill)
- Dedicated drill for each fixture diameter
- F3.5~5.0 drill marking line : bottom line for placing 8.5mm or smaller Fixtures, and top line for 10mm or greater Fixtures
- F5.5 drill marking line : bottom line for placing 6mm or smaller Fixtures, midline for 7mm Fixtures, and top line for 10mm or greater Fixtures
- Drilling up to the lower marking line recommended
- Included in Taper KIT only (not included in 122 Taper KIT)
- F = Fixture




F3.5	F4.0	F4.5	F5.0	F5.5
TCD4C35	TCD4C40	TCD4C45	TCD4C50	TCD4C55

Taper Ultra Drill

09.2013

- Dedicated Taper Drill for Taper Ultra-wide Fixtures of each diameter and length
- Stopper Drill with 1mm margin
- Color coded handle indicating the fixture diameter
- F = Fixture

L	F6.0	F7.0
6	TPD3C 6006	TPD3C 7006
7	TPD3C 6007	TPD3C 7007
8.5	TPD3C 6008	TPD3C 7008
10	TPD3C 6010	TPD3C 7010
11.5	TPD3C 6011	TPD3C 7011
13	TPD3C 6013	TPD3C 7013
Color	Green	Blue




Cortical Drill (Ultra-wide)

01.2009

- Drill used for removing cortical bone from hard bone (for Ultra-wide)
- Dedicated drill for each fixture diameter
- Drilling up to the lower marking line recommended
- F = Fixture


F6.0	F7.0
CD4C60	CD4C70



Parallel Pin (Taper Drill)

- Dedicated Parallel Pin for Taper Drill
- Used for checking the position and direction of bone preparation
- For lower part fixture diameter drill, for upper part Initial Drill
- Color coded according to the fixture diameter (F3.5 : Yellow, F4.0 : Green, F4.5 : Blue, F5.0 : Silver)
- Common component of 122 Taper KIT & Taper KIT

F3.5	F4.0	F4.5	F5.0
TPP3522	TPP4022	TPP4522	TPP5022




Tapered Fixture Tap
(Taper Fixture TSIII, USIII, SSIII SA)

- Dedicated tap for tapered fixture (III type)
- Used for hard bones, forming fixture thread shape
- Torque wrench used after connecting to the engine (25rpm recommended) or a mount extension
- Tapping up to the bottom marking line recommended (F5.0 : Bottom line for placing 7.0mm or smaller Fixtures, and top line for 8.5mm or greater Fixtures)
- F = Fixture

F3.5	F4.0	F4.5	F5.0
OFTS35	OFTS40	OFTS45	OFTS50

※ Refer to surgical instruments for other components (from p142)




Drilling Sequence

Taper Drill

TSIII | SSIII | USIII | KSIII


(Length : 10mm)

Ø 3.0mm




Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Cortical Drill (F3.0)	Ø3.0 Fixture
Soft	☑			Implant Placement
Normal	☑	☑		
Hard	☑	☑	☑	

Ø 3.5mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Cortical Drill (F3.5)	Taper Fixture Tap (F3.5)	Ø3.5 Fixture
Soft	☑	☑				Implant Placement
Normal	☑		☑			
Hard	☑		☑	☑		
Hard (Option)	☑		☑		☑	

Ø 4.0mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.0)	Taper Cortical Drill (F4.0)	Taper Fixture Tap (F4.0)	Ø4.0 Fixture
Soft	☑	☑				Implant Placement
Normal	☑	☑	☑			
Hard	☑	☑	☑	☑		
Hard (Option)	☑	☑	☑		☑	

Taper Cortical Drill marking line: Bottom line for placing 8.5mm or greater Fixtures, and top line for 10mm or greater Fixtures


Recommended placement torque ≤ 40Ncm

TS Fixture placed to a depth 1mm deeper than the bone level for normal bone, and to the bone level for soft bone to maintain fixation stability

For fixture tap used in hard bone, engine (25rpm recommended) is used or Torque Wrench is used after assembling mount extension


(F5.0 Fixture Tap : Bottom line for placing 7.0mm or smaller Fixtures, and top line for 8.5mm or greater Fixtures)

Ø 4.5mm




Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.0)	Taper Drill (F4.5)	Taper Cortical Drill (F4.5)	Taper Fixture Tap (F4.5)	Ø4.5 Fixture
Soft	☑	☑	☑				Implant Placement
Normal	☑	☑		☑			
Hard	☑	☑		☑	☑		
Hard (Option)	☑	☑		☑		☑	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.5)	Taper Drill (F5.0)	Taper Cortical Drill (F5.0)	Taper Fixture Tap (F5.0)	Ø5.0 Fixture
Soft	☑	☑	☑				Implant Placement
Normal	☑	☑	☑	☑			
Hard	☑	☑	☑	☑	☑		
Hard (Option)	☑	☑	☑	☑		☑	

Ø 5.5mm



Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Taper Drill (F4.5)	Taper Drill (F5.5)	Taper Fixture Tap (F5.5)	Ø5.5 Fixture
Soft	☑	☑	☑			Implant Placement
Normal	☑	☑	☑	☑		
Hard	☑	☑	☑	☑	☑	

Drilling Sequence Taper Drill

TSIV | USIV
(Length : 10mm)

Ø 4.0mm

Bone Quality	Twist Drill (Ø2.2)	Taper Drill (F3.5)	Ø4.0 Fixture
D4	☒		Implant Placement
Soft	☒	☒	

Ø 4.5mm

Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Drill (F4.0)	Ø4.5 Fixture
D4		☒			Implant Placement
Soft	☒		☒	☒	

Ø 5.0mm

Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø3.0)	Taper Drill (F3.5)	Taper Drill (F4.5)	Ø5.0 Fixture
D4		☒			Implant Placement
Soft	☒		☒	☒	

Drilling Sequence Taper Drill

TSIII Ultra-wide | SSIII Ultra-wide | USIII Ultra-wide
KSIII Ultra-wide
(Length : 10mm)

Ø 6.0mm

Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Taper Drill (F6.0)	Taper Drill (F6.0)	Cortical Drill (F6.0)	Ø6.0 Fixture
Soft	☒	☒	☒	☒		☒			Implant Placement
Normal	☒	☒	☒	☒	☒		☒		
Hard	☒	☒	☒	☒	☒		☒	☒	

Ø 7.0mm

Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Taper Drill (F6.0)	Taper Drill (F7.0)	Taper Drill (F7.0)	Cortical Drill (F7.0)	Ø7.0 Fixture
Soft	☒	☒	☒	☒		☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒	☒		☒		
Hard	☒	☒	☒	☒	☒	☒		☒	☒	

Recommended placement torque ≤ 40Ncm
TS Fixture placed to a depth 1mm deeper than the bone level for normal bone/hard bone, and to the bone level for soft bone to maintain fixation stability

TSIV Ultra-wide

(Length : 10mm)

Ø 6.0mm

Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.8)	Taper Drill (F6.0)	Ø6.0 Fixture
D4	☒			☒		Implant Placement
Soft	☒	☒	☒		☒	

Ø 7.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.8)	Taper Drill (F7.0)	Ø7.0 Fixture
D4	☒			☒		Implant Placement
Soft	☒	☒	☒		☒	



123 Straight Simple KIT

(O123K) RENEWAL 2020

For **TSII / III** **SSII / III** **USII / III** **KSIII**




123 Straight Simple KIT

Surgical Instruments

123 Twist Drill 03.2012









- Straight Drill to reduce the number of drilling (marking drill)
- Color coded 123 Drill handle indicating the fixture diameter and the main fixture used
- Facilitating drilling depth adjustment by assembling a stopper
- Be sure to use a stopper as it could be difficult to control the depth due to excellent cutting force
- F = Fixture



TL \	D1 / D2 / D3			
	F3.5 (Ø2.2 / 3.0)	F4.0 (Ø3.0 / 3.6)	F4.5 (Ø3.0 / 3.6 / 4.1)	F5.0 (Ø3.0 / 4.1 / 4.6)
34	2D2230FNS	2D3036FNS	2D3041FNS	2D3046FNS
40.4	2D2230FNL	2D3036FNL	2D3041FNL	2D3046FNL
Color	Yellow	Green	Blue	Red

123 Drill Stopper 03.2012

- Number on the stopper indicating the protruding length of the tip when assembled to a drill or instrument
- Color coded by length for easy estimation of the length and relocation of the KIT

L	6.2	7	8	9.5	11	12.5	14	16
								
	ODST05	ODST06	ODST07	ODST08	ODST10	ODST11	ODST13	ODST15


123 Cortical Drill

10.2011


- Drill used for removing cortical bone from hard bone
- Drilling up to the bottom marking line recommended
- II type marking line : for hard bone
- III type marking line : bottom line for normal bone, and top line for hard bone
- IV type marking line : for normal bone
- Color coded handle indicating the fixture diameter and the main fixture used
- F = Fixture

Type	F3.5	F4.0	F4.5	F5.0
II	O2CD 35	O2CD 40	O2CD 45	O2CD 50
III	O3CD 35	O3CD 40	O3CD 45	O3CD 50
Color	Yellow	Green	Blue	Red


※ Refer to surgical instruments for other components (from p142)



123 Straight KIT (O123FK) RENEWAL 2020




123 Straight Full KIT (O123STFK) 04.2018



123 Straight KIT Surgical Instruments

123 Twist Drill (Stopper Drill) 06.2013

- Straight Drill to reduce the number of drilling (with stopper)
- Color coded 123 Drill handle indicating the fixture diameter and the main fixture used
- F = Fixture




L	TL	D1 / D2 / D3			
		F3.5 (Ø2.2 / 3.0)	F4.0 (Ø3.0 / 3.6)	F4.5 (Ø3.0 / 3.6 / 4.1)	F5.0 (Ø3.0 / 4.1 / 4.6)
6	30.5	2D2230 06LC	2D3036 06LC	2D3041 06LC	2D3046 06LC
7	31.5	2D2230 07LC	2D3036 07LC	2D3041 07LC	2D3046 07LC
8.5	33	2D2230 08LC	2D3036 08LC	2D3041 08LC	2D3046 08LC
10	34.5	2D2230 10LC	2D3036 10LC	2D3041 10LC	2D3046 10LC
11.5	34.5	2D2230 11LC	2D3036 11LC	2D3041 11LC	2D3046 11LC
13	36	2D2230 13LC	2D3036 13LC	2D3041 13LC	2D3046 13LC
15	38	2D2230 15LC	2D3036 15LC	2D3041 15LC	2D3046 15LC
Color		Yellow	Green	Blue	Red

123 Cortical Drill 10.2011

- Drill used for removing cortical bone from hard bone
- Drilling up to the bottom marking line recommended
- II type marking line : for hard bone
- III type marking line : bottom line for normal bone, and top line for hard bone
- IV type marking line : for normal bone
- Color coded handle indicating the fixture diameter and the main fixture used
- F = Fixture


Type	F3.5	F4.0	F4.5	F5.0
II	O2CD 35	O2CD 40	O2CD 45	O2CD 50
III	O3CD 35	O3CD 40	O3CD 45	O3CD 50
Color	Yellow	Green	Blue	Red



Parallel Pin (123 Drill) 03.2012


- Dedicated Parallel Pin for 123 Twist Drill
- Used for checking the position and direction of bone preparation
- Bottom part for Initial Drill, and top part for F3.5(Ø2.2/3.0) drill

D	Ø4.0	Ø 5.0
	OPLP400	OPLP500



123 Ultra Twist Drill

- 2-stage drill with both Pilot and Twist Drill functions
- Straight Drill to reduce the number of drilling (with stopper)
- Dedicated drill used for F7.0 Fixtures in soft bone
- F = Fixture




L	TL	F5.5 (Ø4.6 / 5.2)	F6.0 (Ø4.6 / 5.5)	F7.0Soft (Ø5.5 / 6.2)	F7.0 (Ø5.5 / 6.5)
6	30.5	3D4652 06T	3D4655 06T	-	3D5565 06T
7	31.5	3D4652 07T	3D4655 07T	-	3D5565 07T
8.5	33.5	3D4652 08T	3D4655 08T	-	3D5565 08T
10	34.5	3D4652 10T	3D4655 10T	-	3D5565 10T
11.5	34.5	3D4652 11T	3D4655 11T	-	3D5565 11T
13	36.0	3D4652 13T	3D4655 13T	3D5562 13T	3D5565 13T
Color		Yellow	Green	Blue	Red

Trial Pin (Ultra-wide) 01.2009

- Checking the width and depth of a fresh extraction socket or failed implant socket
- Checking the drilling after using a Direct Drill as the final drill
- Used as a Parallel Pin

D	Ø 5.2	Ø 5.5	Ø 6.2	Ø 6.5
	UWFTP52	UWFTP55	UWFTP62	UWFTP65



※ Refer to surgical instruments for other components (from p142)


Drilling Sequence II Type 123 Twist Drill

TSII | SSII | USII

(Length : 10mm)




Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Cortical Drill	Ø3.5 Fixture
Soft	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		Implant Placement
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.1)	Twist Drill (Ø3.0/4.6)	Cortical Drill	Ø5.0 Fixture
Soft	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>			Implant Placement
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Cortical Drill	Ø4.0 Fixture
Soft	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			Implant Placement
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Twist Drill (Ø3.0/4.1)	Cortical Drill	Ø4.5 Fixture
Soft	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>			Implant Placement
Normal	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>		
Hard	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	


Drilling Sequence III Type 123 Twist Drill

TSIII | SSIII | USIII | KSIII

(Length : 10mm)



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Cortical Drill (F3.5) Bottom line	Cortical Drill (F3.5) Top line	Ø3.5 Fixture
Soft	☒	☒			Implant Placement
Normal	☒	☒	☒		
Hard	☒	☒		☒	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.1)	Twist Drill (Ø3.0/4.6)	Cortical Drill (F5.0) Bottom line	Cortical Drill (F5.0) Top line	Ø5.0 Fixture
Soft	☒		☒				Implant Placement
Normal	☒	☒		☒	☒		
Hard	☒	☒		☒		☒	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Cortical Drill (F4.0) Bottom line	Cortical Drill (F4.0) Top line	Ø4.0 Fixture
Soft	☒	☒				Implant Placement
Normal	☒	☒	☒	☒		
Hard	☒	☒	☒		☒	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Cortical Drill (F5.5) Bottom line	Cortical Drill (F5.5) Top line	Ø5.5 Fixture
Soft	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒		
Hard	☒	☒	☒		☒	



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/3.6)	Twist Drill (Ø3.0/4.1)	Cortical Drill (F4.5) Bottom line	Cortical Drill (F5.0) Top line	Ø4.5 Fixture
Soft	☒		☒				Implant Placement
Normal	☒	☒		☒	☒		
Hard	☒	☒		☒		☒	

Drilling Sequence Ultra-wide 123 Twist Drill

TSII Ultra-wide | SSII Ultra-wide | USII Ultra-wide

(Length : 10mm)

Ø 6.0mm

Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.2)	Twist Drill (Ø4.6/5.5)	Cortical Drill (F6.0)	Ø 6.0 Fixture
Soft	☒	☒	☒	☒			
Normal	☒	☒	☒		☒		Implant Placement
Hard	☒	☒	☒		☒	☒	

Ø 7.0mm

Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.5)	Twist Drill (Ø4.6/5.5) (F7.0 Soft)	Twist Drill (Ø5.5/6.5)	Cortical Drill (F7.0)	Ø 7.0 Fixture
Soft	☒	☒	☒	☒	☒			
Normal	☒	☒	☒	☒			☒	Implant Placement
Hard	☒	☒	☒	☒		☒	☒	

Drilling Sequence Ultra-wide 123 Twist Drill

TSIII Ultra-wide | SSIII Ultra-wide | USIII Ultra-wide

KSIII Ultra-wide

(Length : 10mm)

Ø 6.0mm



Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.2)	Twist Drill (Ø4.6/5.5)	Cortical Drill (F6.0)	Ø 6.0 Fixture
Soft	☒	☒	☒	☒(6mm)			
Normal	☒	☒	☒		☒(6mm)		Implant Placement
Hard	☒	☒	☒		☒(6mm)	☒	

Ø 7.0mm

Bone Quality	Guide Drill	Twist Drill (Ø2.2/3.0)	Twist Drill (Ø3.0/4.6)	Twist Drill (Ø4.6/5.5)	Twist Drill (Ø5.5/6.2) (F7.0 Soft)	Twist Drill (Ø5.5/6.5)	Cortical Drill (F7.0)	Ø 7.0 Fixture
Soft	☒	☒	☒	☒	☒(6mm)			
Normal	☒	☒	☒	☒			☒(6mm)	Implant Placement
Hard	☒	☒	☒	☒		☒(6mm)	☒	


New Hanaro KIT (HKA2) 03.2013

For **TSII / III** **SSII / III** **USII / III** **KSIII**



Twist Drill (Stopper Drill) 12.2012


- Long stopper (6mm)
- Enabling a procedure without drill
- Color coded stopper indicating the drill length



L	TL	D	Ø2.2	Ø3.0	Ø3.3	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
	Y-Dim		0.6	0.9	1.0	1.0	1.0	1.0	1.0	1.0
6	30.5	2D22 06LC	3D30 06LC	-	-	3D38 06LC	-	-	-	-
7	31.5	2D22 07LC01	3D30 07LC01	-	-	3D38 07LC01	-	-	-	-
8.5	33	2D22 08LC01	3D30 08LC01	-	-	3D38 08LC01	-	-	-	-
10	34.5	2D22 10LC01	3D30 10LC01	-	-	3D38 10LC01	-	-	-	-
11.5	34.5	2D22 11LC01	3D30 11LC01	3D33 11LC01	3D36 11LC01	3D38 11LC01	3D41 11LC01	3D43 11LC01	3D46 11LC01	-
13	36	2D22 13LC01	3D30 13LC01	-	-	3D38 13LC01	-	-	-	-

Twist Drill (Non-Stopper Drill) 01.2009

- Used for limited Stopper Drill access into the oral cavity
- Refer to the Non-stopper Drill image for marking drill marking line sizes for Short/Long types




TL	D	Ø1.5	Ø2.0	Ø2.2	Ø2.7	Ø3.0	Ø3.3
33		2D15 18FNLC	2D20 18FNLC	2D22 18FNLC	3D27 18FNLC	3D30 18FNLC	3D33 18FNLC
41		-	-	2D22 15FNLC01	3D27 15FNLC01	3D30 15FNLC01	3D33 15FNLC01

TL	D	Ø3.6	Ø3.8	Ø4.1	Ø4.3	Ø4.6
33		3D36 18FNLC	3D38 18FNLC	3D41 18FNLC	3D43 18FNLC	3D46 18FNLC
41		3D36 15FNLC01	3D38 15FNLC01	3D41 15FNLC01	3D43 15FNLC01	3D46 15FNLC01

Long Shank Pilot Drill 01.2009

- Used for adjusting the drilling hole path
- Previous drilling path maintained for the next drill

D1 / D2	Ø2.0 / 2.7	Ø2.0 / 3.0	Ø3.0 / 3.8	Ø3.0 / 4.1
	APD270C	APD300C	APD380C	APD410C



Cortical Drill 2 (TSII, SSII SA) 01.2009

- Drill used for removing cortical bone from hard bone (for II type)
- Dedicated drills available for each fixture diameter
- Drilling up to the bottom marking line recommended
- F = Fixture

F3.5	F4.0	F4.5	F5.0
CD2C35	CD2C40	CD2C45	CD2C50



Cortical Drill 3 (Taper Fixture TSIII, SSIII, USIII, KSIII) 08.2014

- Drill used for expanding the cortical bone after using the Straight Drill
- Used after forming the final drill hole in normal or harder bone
- Dedicated drills available for each fixture diameter
- Bottom marking line for normal bone, and top marking line for hard bone
- Drilling up to the lower marking line recommended

F3.0	F3.5	F4.0	F4.5	F5.0	F5.5
CD4C30	CD4C35	CD4C40	CD4C45	CD4C50	CD4C55



Countersink (USIII, USII SA, USIII SA Wide PS, Wide) 01.2009

- Dedicated drill for expanding the placement hole opening for US Fixtures
- wide PS and wide of USIII, USII SA, and USIII SA
- Recommended speed : 300rpm

USSCS45W



Straight Fixture Tap

(TSII, USII, SSII SA)

02.2016

- Dedicated tap for Straight Fixtures (II type)
- Used for hard bones, forming fixture thread shape
- Torque wrench used after connecting to the engine (25rpm recommended)
or a mount extension
- Tapping up to the bottom marking line recommended
- F = Fixture

F3.5	F4.0	F4.5	F5.0
O2FTS35	O2FTS40	O2FTS45	O2FTS50




Parallel Pin

01.2013

- Used for checking the position and direction of bone preparation

D	Ø4.0	Ø5.0	Ø6.0	Full Set
	APP400	APP500	APP600	APPS

※ Refer to surgical instruments for other components (from p142)




For Ultra-wide

Lower panel components

Open Wrench
SPOW




Ratchet Wrench
RCWC



Direct Drill 01.2009


- 2-stage drill with both Pilot and Twist Drill functions
- Enabling final drilling without pilot drilling
- Increased initial fixation stability in a fresh extraction socket with reduced dead space in apex



D1 / D2	Ø4.6 / 5.2	Ø4.6 / 5.5	Ø5.5 / 6.2	Ø5.5 / 6.5
	3D5213FNLC	3D5513FNLC	3D6213FNLC	3D6513FNLC

Cortical Drill (Ultra-wide) 01.2009


- Drill used for removing cortical bone from hard bone (for Ultra-wide)
- Dedicated drills available for each fixture diameter
- Drilling up to the lower marking line recommended
- F = Fixture



	F6.0	F7.0
	CD4C60	CD4C70

Trial Pin (Ultra-wide) 01.2009


- Checking the width and depth of a fresh extraction socket or failed implant socket
- Checking the drilling after using a Direct Drill as the final drill
- Used as a Parallel Pin



D	Ø5.2	Ø5.5	Ø6.2	Ø6.5
	UWFTP52	UWFTP55	UWFTP62	UWFTP65


※ Refer to surgical instruments for other components (from p142)

Ø 3.5mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Cortical Drill (F3.5)	Straight Fixture Tap (F3.5)	Ø3.5 Fixture
Soft	☒	☒					Implant Placement
Normal	☒		☒	☒			
Hard	☒		☒	☒	☒		
Hard (Option)	☒		☒	☒		☒	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.3)	Twist Drill (Ø4.6)	Cortical Drill (F5.0)	Straight Fixture Tap (F5.0)	Ø5.0 Fixture
Soft	☒	☒	☒	☒	☒	☒				Implant Placement
Normal	☒	☒	☒	☒	☒		☒			
Hard	☒	☒	☒	☒	☒		☒	☒		
Hard (Option)	☒	☒	☒	☒	☒		☒		☒	

Ø 4.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.3)	Twist Drill (Ø3.8)	Cortical Drill (F4.0)	Straight Fixture Tap (F4.0)	Ø4.0 Fixture
Soft	☒	☒	☒	☒				Implant Placement
Normal	☒	☒	☒		☒			
Hard	☒	☒	☒		☒	☒		
Hard (Option)	☒	☒	☒		☒		☒	

Ø 4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.1)	Cortical Drill (F4.5)	Straight Fixture Tap (F4.5)	Ø4.5 Fixture
Soft	☒	☒	☒	☒	☒				Implant Placement
Normal	☒	☒	☒	☒	☒	☒			
Hard	☒	☒	☒	☒	☒	☒	☒		
Hard (Option)	☒	☒	☒	☒	☒	☒		☒	

Recommended placement torque

≤ 40Ncm

TS Fixture placed to a depth

1mm deeper than the bone level for normal bone/hard bone, and to the bone level for soft bone to maintain fixation stability

For fixture tap used in hard bone, engine (25rpm recommended) is used or Torque Wrench is used after assembling mount extension

Drilling Sequence III Type Straight Drill

TSIII | SSIII | USIII | KSIII

(Length : 10mm)

Ø 3.0mm



Bone Quality	Twist Drill (Ø2.2)	Twist Drill (Ø2.7)	Cortical Drill 2 (F3.0)	Ø3.0 Fixture
Soft	☒			Implant Placement
Normal	☒	☒		
Hard	☒	☒	☒	

Ø 3.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø2.2)	Cortical Drill 3 (F3.5)	Cortical Drill 3 (F3.5)	Ø3.5 Fixture
Soft	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒		
Hard	☒	☒	☒		☒	

Ø 4.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.3)	Cortical Drill 3 (F4.0)	Cortical Drill 3 (F4.0)	Ø4.0 Fixture
Soft	☒	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒		
Hard	☒	☒	☒	☒		☒	

Recommended placement torque ≤ 40Ncm

TS Fixture placed to a depth 1mm deeper than the bone level for normal bone/hard bone, and to the bone level for soft bone to maintain fixation stability

Ø 4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Cortical Drill 3 (F4.5)	Cortical Drill 3 (F4.5)	Ø4.5 Fixture
Soft	☒	☒	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒	☒		
Hard	☒	☒	☒	☒	☒		☒	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.3)	Cortical Drill 3 (F5.0)	Cortical Drill 3 (F5.0)	Ø5.0 Fixture
Soft	☒	☒	☒	☒	☒				Implant Placement
Normal	☒	☒	☒	☒	☒	☒	☒		
Hard	☒	☒	☒	☒	☒	☒		☒	

Ø 5.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Cortical Drill 3 (F5.5)	Cortical Drill 3 (F5.5)	Ø5.5 Fixture
Soft	☒	☒	☒	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒	☒	☒		
Hard	☒	☒	☒	☒	☒	☒		☒	

Drilling Sequence IV Type Straight Drill

TSIV | USIV
(Length : 10mm)

Ø 4.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Twist Drill (Ø3.0 Half)	Ø4.0 Fixture
D4	☒				Implant Placement
Soft	☒	☒	☒	☒	

Ø 4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.1 Half)	Ø4.5 Fixture
D4			☒				Implant Placement
Soft	☒	☒	☒	☒	☒	☒	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6 Half)	Ø5.0 Fixture
D4			☒				Implant Placement
Soft	☒	☒	☒	☒	☒	☒	

Drilling Sequence Ultra-wide Straight Drill

TSII Ultra-wide | SSII Ultra-wide | USII Ultra-wide
(Length : 10mm)

Ø 6.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Direct Drill (Ø5.2)	Direct Drill (Ø5.5)	Cortical Drill (F6.0)	Ø6.0 Fixture
Soft	☒	☒	☒	☒	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒	☒		☒		
Hard	☒	☒	☒	☒	☒	☒		☒	☒	

Ø 7.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Direct Drill (Ø5.5)	Direct Drill (Ø6.2)	Direct Drill (Ø6.5)	Cortical Drill (F7.0)	Ø7.0 Fixture
Soft	☒	☒	☒	☒	☒	☒	☒	☒			Implant Placement
Normal	☒	☒	☒	☒	☒	☒	☒		☒		
Hard	☒	☒	☒	☒	☒	☒	☒		☒	☒	

Recommended placement torque ≤ 40Ncm
TSIV/USIV Fixtures are dedicated implants for maxillary sinus or soft bone, not guiding normal or harder bones
Reducing the speed to 15rpm or lower recommended for placement as the placement speed is too fast for TSIV/USIV Fixtures due to large thread pitch

Drilling Sequence

Ultra-wide

Straight Drill

TSII Ultra-wide

SSII Ultra-wide

USII Ultra-wide

(Length : 10mm)

Ø 6.0mm
Drilling sequence with trephine in the healed mature bone



10mm

Bone Quality	Trephine Drill (Ø4.2/5.0)	Direct Drill (Ø4.6/5.2)	Direct Drill (Ø4.6/5.5)	Cortical Drill (F6.0)	Ø6.0 Fixture
Soft	☒	☒			
Normal	☒		☒		Implant Placement
Hard	☒		☒	☒	

Immediate placement at the extraction socket



10mm

Trephine Drill (Ø4.2/5.0) Direct Drill (Ø4.6/5.0) Trial Pin (Ø5.5) 6.0 Ultra-wide Fixture

Immediate replacement of the failed implant



10mm

Ø4.0 Failed Fixture Trephine Drill (Ø4.2/Ø5.0) Trial Pin (Ø5.5) Direct Drill (Ø4.6/Ø5.5) 6.0 Ultra-wide Fixture

Drilling Sequence

Ultra-wide

Straight Drill

TSIII Ultra-wide

SSIII Ultra-wide

USIII Ultra-wide

KSIII Ultra-wide

(Length : 10mm)

Ø 6.0mm



10mm

Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Direct Drill (Ø5.2)	Direct Drill (Ø5.5)	Cortical Drill (F6.0)	Ø6.0 Fixture
Soft	☒	☒		☒	☒	☒	☒			
Normal	☒	☒		☒	☒	☒		☒		Implant Placement
Hard	☒	☒		☒	☒	☒		☒	☒	

Ø 7.0mm



10mm

Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.0)	Twist Drill (Ø4.6)	Direct Drill (Ø5.5)	Direct Drill (Ø6.2)	Direct Drill (Ø6.5)	Cortical Drill (F7.0)	Ø7.0 Fixture
Soft	☒	☒	☒	☒	☒	☒	☒	☒			
Normal	☒	☒	☒	☒	☒	☒	☒		☒		Implant Placement
Hard	☒	☒	☒	☒	☒	☒	☒		☒		

Recommended placement torque ≤ 40Ncm
TS Fixture placed to a depth 1mm deeper than the bone level for normal bone/hard bone, and to the bone level for soft bone to maintain fixation stability

Drilling Sequence

Ultra-wide


Straight Drill

TSIV Ultra-wide

USIV Ultra-wide

(Length : 10mm)

Ø 6.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Direct Drill (Ø5.2)	Ø6.0 Fixture
D4	☒	☒			☒			Implant Placement
Soft	☒	☒	☒		☒		☒	

Ø 7.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill (Ø2.0/3.0)	Twist Drill (Ø3.0)	Pilot Drill (Ø3.0/3.8)	Twist Drill (Ø3.8)	Twist Drill (Ø4.6)	Direct Drill (Ø5.5)	Direct Drill (Ø6.2)	Ø7.0 Fixture
D4	☒	☒			☒	☒			Implant Placement
Soft	☒	☒		☒	☒	☒	☒	☒	

OSSTEM[®]
IMPLANT


Recommended placement torque ≤ 40Ncm

For TSIII SSIII USIII KSIII

Lower panel components


Torque Wrench
TW30B

Depth Gauge
OSDG



485 Drill

- Drill for placing short implants in alveolar bone lacking in vertical height
- Ø2.2 drill : Straight Drill
- Top blade of other drill in the shape of CAS Drill, and the side blade in the shape of Taper Drill
- Stopper Drill with 1mm margin
- Recommended speed : 800~1,200rpm



L	Type	Ø2.2	Pilot	F4.0	F4.5	F5.0	F5.5
4.0		O485D 2204	O485D 3504	O485D 4004	O485D 4504	O485D 5004	O485D 5504
5.0		O485D 2205	-	O485D 4005	O485D 4505	O485D 5005	O485D 5505
6.0		O485D 2206	O485D 3506	O485D 4006	O485D 4506	O485D 5006	O485D 5506
7.0		O485D 2207	-	O485D 4007	O485D 4507	O485D 5007	O485D 5507
8.5		O485D 2208	-	O485D 4008	O485D 4508	O485D 5008	O485D 5508


※ Refer to surgical instruments for other components (from p142)

Drilling Sequence 485 Drill

TSIII | SSIII | USIII | KSIII


(Length : 7mm)

Ø 4.0mm




Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.0)	485 Drill (F4.5)	Ø4.0 Fixture
Normal	☒	☒	☒		Implant Placement
Hard	☒	☒		☒	

Ø 4.5mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F4.5)	485 Drill (F5.0)	Ø4.5 Fixture
Normal	☒	☒	☒		Implant Placement
Hard	☒	☒		☒	

Ø 5.0mm



Bone Quality	Twist Drill (Ø2.2)	Pilot Drill	485 Drill (F5.0)	485 Drill (F5.5)	Ø5.0 Fixture
Normal	☒	☒	☒		Implant Placement
Hard	☒	☒		☒	


OSSTEM[®]
IMPLANT

Surgical Instruments

123 Guide Drill

- Drill for forming a hole to facilitate initial drilling
- Facilitating drilling depth adjustment by assembling a stopper
- Included in 122 Taper KIT only (not included in Taper KIT)

D	Ø 2.0
OGD2027L	



Drill Extension

- Extending the length of a drill or other hand piece tool (drill extended by 16.9mm)
- Risk of bending or fracture upon exerting excessive force on inadequate assembly
- Common component of Taper KIT and Straight KIT


L (Extention)	16.9
ODE	



Lance Drill (Guide Drill)

- Forming a hole to facilitate initial drilling
- Bone density determined through drilling
- Included in Taper KIT only (not included in 122 Taper KIT)

L	Short	Long
AGDSC		AGDLC




TS NoMount Driver

05.2012

- Driver directly connected to the fixture upon placing with a surgical hand piece
- C = Connection


L \ C	Mini	Regular
Short	TSNMDMS	TSNMDS
Long	TSNMDML	TSNMDRL
Ex.Long	TSNMDME	TSNMDRE



Sidecut Drill

- Drill to remove the side parts with the cutting edge of the body
- Used to remove the ridge of a fresh extraction socket
- Facilitating site preparation of a fresh extraction socket
- Included in Taper KIT only (not included in 122 Taper KIT)

L \ D1 / D2	Ø 1.5 / 2.0	Ø 2.0 / 2.5	Ø 3.0 / 3.5
13	OSLM DS	OSLMD 20S	-
16.5	-	-	OSLMD 30L
20	OSLM DL	OSLMD 20L	-




SS NoMount Driver

09.2014

- Driver directly connected to the fixture upon placing with a surgical hand piece
- C = Connection

L \ C	Regular / Wide
Short	SSNMDS
Long	SSNMDL
Ex.Long	SSNMDE




Surgical Instruments

US NoMount Driver 12.2009

- Driver directly connected to the fixture upon placing with a surgical hand piece
- C = Connection


L \ C	Mini	Regular	Wide
Short	USNMD35MS	USNMD41RS	USNMD51WS
Long	USNMD35ML	USNMD41RL	USNMD51WL



SS NoMount Torque Driver 01.2009

- Driver directly connected to the fixture upon placing with a wrench
- Be sure to use it after confirming an adequate assembly (Risk of fracture even at low torque when inadequately assembled)
- Note that it cannot be removed in case of fracture
- C = Connection


L \ C	Regular / Wide
Short	SSNMT39S
Long	SSNMT39L



KS NoMount Driver 10.2019

- Driver directly connected to the fixture upon placing with a surgical hand piece
- Ø3.5 Fixtures assembled below the bottom marking; and Ø4.0, Ø4.5, Ø5.0, Ø6.0 and Ø7.0 Fixtures assembled above the bottom marking
- Distance between laser markings and laser marking are divided into 0.5mm
- C = Connection


L \ C	Regular
Short	KSNMDS
Long	KSNMDL



TS Fixture Driver 11.2014

- Used by assembling directly to the fixture for final placement depth adjustment or removal
- C = Connection


L \ C	Mini	Regular
Short	GSMFDS	GSRFDS
Long	GSMFDL	GSRFDL
Ex.Long	GSMFDE	GSRFDE



TS NoMount Torque Driver 12.2009

- Driver directly connected to the fixture upon placing with a wrench
- Be sure to use it after confirming an adequate assembly (Risk of fracture even at low torque when inadequately assembled)
- Note that it cannot be removed in case of fracture
- C = Connection


L \ C	Mini	Regular
Short	GSNMT32S	GSNMT35S
Long	GSNMT32L	GSNMT35L
Ex.Long	GSNMT32E	GSNMT35E



SS Fixture Driver 12.2014

- Used by assembling directly to the fixture for final placement depth adjustment or removal
- C = Connection

L \ C	Regular / Wide
Short	SSRFDS
Long	SSRFDL
Ex.Long	SSRFDE



Surgical Instruments

US Fixture Driver 01.2009

- Used by assembling directly to the fixture for final placement depth adjustment or removal
- C = Connection


C	Mini	Regular	Wide
	USMFDL	USRFDL	USWFDL



Simple Mount Driver 01.2009

- Used by assembling to the simple mount for fixture placement

L	
Short	ASMDS
Long	ASMDL



KS Fixture Driver 10.2019

- Used by assembling directly to the fixture for final placement depth adjustment or removal
- C= Connection


L	C	Regular
Short		KSFDS
Long		KSFDL



Simple Mount Extension 12.2014

- Used by connecting to a wrench for extending the simple mount length or applying torque manually

L	
Short	ASMES
Long	ASMEL



Torque Extension 12.2013

- Extending the length of the instrument used by connecting to a wrench (10mm extension)

OTE



Simple Open Wrench 01.2009

- Used for removing a simple mount from weak bone
- Easy placement into the oral cavity with 30°

ASOW



Surgical Instruments

Removal Tool (Fixture Mount) 01.2009

- Used after removing mount screw in case of jamming between the fixture and mount
- Used by assembling to driver handle and Torque Wrench
- Removing mount by rotating FWD after inserting vertically
- App = Application



App	Mini (TS,US)	Regular (TS,SS,US)	Wide (SS)	Wide (US)
	ERFM	HRFR		ERFW

Depth Gauge

- Used for measuring the drilling depth (7-15mm) or as an open wrench
- Common component of 122 Taper & Taper KIT




	OSDG
--	------

Positioning Guide 01.2009

- Instrument to facilitate drilling interval setting for fixture placement
- Placed into the hole for use after initial drilling
- Packing unit : each component or the set

W / L	2.5 / 21.5	6.0 / 17.5	11 / 17.5
	APG201	APG202	APG203



Tissue Height Gauge (TS) 01.2009

- Instrument to measure the gingival height by assembling to the fixture connection for top G/H selection in TS implant placement

	HGTSHG
--	--------



Ratchet Wrench 01.2009

- Dedicated wrench for anti-reverse procedure
- Excessive torque exertion may result in internal damage to bone or fixture


	CITQW-1185A
--	-------------



L-Wrench 10.2013

- 1.2 hex driver for overcoming narrow spacing
- Torque indication : 5~8Ncm torque at the level when the wrench appears to be bent a little (within 10°)

	LWC
--	-----



Torque Wrench (Spring Type) 06.2012

- Wrench to apply a constant torque (10/20/30Ncm) to screws and abutments
- When the set torque is applied, the neck of the Torque Wrench is bent for indication
- If a continuous force is applied while the neck is bent, excessive torque is applied, resulting in screw fracture

	TW30
--	------



Torque Wrench (Bar Type) 05.2012

- Used for adjusting the implant placement position and tightening screws and abutments
- Applying torque according to the line marked with the torque value to be applied by pulling the bar

	TW30B
--	-------



Surgical Instruments

Torque Wrench Set 11.2015


- A set of a two-way Torque Wrench and a Torque Connector
- Applying forward/reverse torque by rotating the Torque Wrench handle without removing the connector
- Compatible with osstem machine driver connector
- Applying torque according to the line marked with the torque value to be applied by pulling the bar
- Packing unit : changeable Torque Wrench + Torque Connector



MX30

Tissue Punch 09.2011

- Instrument used for flapless surgery
- Marked at 2mm intervals for measuring gingival height
- Packing unit : Tissue Punch + Guide pin
- ※ Using a Tissue Punch with a smaller diameter than the Healing Abutment recommended




D	Ø3.3	Ø3.8	Ø4.3	Ø4.8	Ø5.3
	OSTP33	OSTP38	OSTP43	OSTP48	OSTP53
TS	Ø 4.0/4.5	Ø 4.5/5.0	Ø 5.0	Ø 6.0	Ø 6.0
SS	-	Ø 4.8	-	Ø 6.0	Ø 6.0
US	Ø 4.0	Ø 5.0	Ø 5.0	Ø 6.0	Ø 6.0

For application Healing Abutment

Bone Profiler (TS) 01.2009

- Used for removing bone around the fixture for the 1st and 2nd stage surgery
- Used by connecting a guide screw to the fixture and removing bone to compensate for the shape of the Healing Abutment
- Guide Screw protecting the morse taper entrance of the fixture
- Packing unit : Bone Profiler + Guide Screw
- C = Connection



C \ D (Healing Abutment)	Ø4.5	Ø5.5	Ø6.5 / 7.5
Mini / Regular	GSBP45	GSBP55	GSBP75
	Mini + Regular Guide Screw	Mini + Regular Guide Screw	Regular Guide Screw

Bone Profiler (US) 01.2009

- Used for removing bone formed around the cover screw for the secondary procedure
- Used by connecting a guide screw to the fixture to compensate for the shape of the Healing Abutment
- Guide Screw protecting the fixture hex
- Packing unit : Bone Profiler + Guide Screw
- P = Platform




D \ P	Mini	Regular	Wide	T-type
Ø4.0	ABPM 400C	-	-	-
Ø5.0	ABPM 500C	ABPR 500C	-	-
Ø6.0	-	ABPR 600C	ABPW 600C	TBPW 600C
Ø7.0	-	-	ABPW 700C	-

Surgical Instruments

Trephine Drill 01.2009

- Used for collecting bone or for removing damaged or failed fixtures
- Used for removing septal bone
- Used as an Initial Drill for ultra-fixtue placement



L \ D (Inner / Outer)	3.7 / 4.5	4.2 / 5.0	4.7 / 5.5	5.2 / 6.0	5.7 / 6.5	6.2 / 7.0
Short	TD37S	TD42S	TD47S	TD52S	TD57S	TD62S
Long	TD37	TD42	TD47	TD52	TD57	TD62


Anterior Hand Driver (Implant) 12.2014

- Instrument for manual placement in anterior region
- Used by connecting to a NoMount Torque Driver or Fixture Driver
- Excessive torque may result in fracture of the fixture or driver



Machine Driver Handle 12.2013

- Enabling hand rotation by connecting to any surgical instrument for engine




Torque Handle 11.2015

- Manual instrument used by connecting to the contra-angle hand piece (1:1 gear ratio for hand piece)
- Used for tightening screws such as Healing Abutment, Cover Screw, Abutment Screw and Orthodontic Screw (used for temporary tightening of Abutment Screw, which requires final tightening with a Torque Wrench)
- Excessive torque may result in fracture or malfunction of the hand piece



Bone Mill 01.2009

- Forming particulate bone with collected autogenous bone



Prosthetic Simple KIT (OPSK) 02.2017

Prosthetic KIT (OPK) 05.2018


Top panel components

Torque Wrench
TW30B




154

OSSTEM KIT




155

OSSTEM KIT



Hand Driver


- Manual driver
- Tip holding feature (except internal hex type)
- Internal hex type length : 11



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Ex.Short (8)	-	AHD09MSH	AHD12MSH	-	-
Short (13)	ASD 05SH	AHD09SH	AHD12SH	IHD20H	IHD27H
Middle (15)	-	-	AHD12MH	-	-
Long (18)	ASD 05LH	AHD09LH	AHD12LH	-	-
Ex.Long (25)	-	-	AHD12EH	-	-

Torque Driver


- Driver for Torque Wrench assembly
- Tip holding feature
- Use the recommended torque (excessive torque may result in fracture)
- Risk of fracture even at low torque when inadequately assembled
- Exerting torque with the driver straight up (with no tilting)
- Be sure to replace any bent tips due to extended use or excessive torque



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Ex.Short (8)	-	-	TRHD 12MS	-	-
Short (13)	TRSD 05S	TRHD 09S	TRHD 12S	TIHD20S	-
Middle (15)	-	-	TRHD 12M	-	-
Long (20)	TRSD 05L	TRHD 09L	TRHD 12L	TIHD20L	TIHD27
Ex.Long (25)	TRSD 05E	-	TRHD 12E	-	-

Machine Screw Driver

- Driver for engine
- Tip holding Tip holding feature (except internal hex type)
- Internal hex type length : 8



L \ Type	0.5 Slot	0.9 Hex	1.2 Hex	2.0 Int.Hex	2.7 Int.Hex
Osstem Torque (5)	-	-	OTH12S	-	-
Short (5.6)	AMSD 05S	AMSD 09S	AMSD 12S	-	-
Long (11.6)	AMSD 05L	AMSD 09L	AMSD 12L	EIHD20	EIHD27
Ex.Long (17.6)	-	-	AMSD 12E	-	-


Application

Product applied to a driver
(Common use for hand, Machine Screw,
and Torque Driver)

Cover Screw (US Mini)	Healing Abutment, UCLA, Cemented Abutment Screw, Mount Screw	Esthetic Abutment Screw Regular, Esthetic-low Abutment Screw, Standard	Wide Esthetic-low Abutment Screw
--------------------------	--	---	-------------------------------------

Angled Torque Driver 02.2017


- Driver for Torque Wrench assembly
- No holding feature
- Recommended tightening torque : 30Ncm (excessive torque may result in fracture)
- Do not remove the tube preventing debris upon fracture
- Recommended use cycle : 10 times
- Set : 3ea



L \ Type	1.2 Hex	1.2 Hex (Set)
Short (13)	ATD12S	ATD12S3S
Long (20)	ATD12L	ATD12L3S


Repair Torque Driver

- Handle diameter reduced compared to Torque Driver (Ø2.1 → Ø1.6)
- Minimizing crown hole diameter for prosthesis repair or SCRP procedure



Solid Abutment Driver

- Dedicated driver for solid abutment
- Applying torque after inserting the groove of the solid abutment to the part with a triangular marking
- Recommended tightening torque : 30Ncm



Regular


L \ Type	Square	Round
Short (6)	SDSS	SDRS
Long (12)	SDSL	SDRL

Wide

L \ Type	Square
Short (10)	SD60S

Excellent Solid Abutment Driver

- Dedicated driver for excellent solid abutment
- Applying torque after inserting the groove of the excellent solid abutment to the part with a triangular marking
- Recommended tightening torque : 30Ncm



Regular

L \ Type	Square	Round
Short (6)	ESDSS	ESDRS
Long (12)	ESDSL	ESDRL


Wide

L \ Type	Square
Short (10)	ESD60S

O-ring Abutment Driver

- Dedicated driver for O-ring Abutment

	Small	Normal
	STAOB	AORB



Rigid Outer Driver

- Dedicated driver for Rigid Abutment
- Recommended tightening torque : 30Ncm


L \ D (Abutment)	Ø4.0	Ø4.5	Ø5.0	Ø6.0
Short (16.5)	ORDMS	ORD45S	ORDRS	ORDWS
Long (21.5)	ORDML	ORD45L	ORDRL	ORDWL



Octa Abutment Driver

- Dedicated driver for Octa Abutment
- Recommended tightening torque : 30Ncm

L \ Type	Square	Round
Short	ODSS	ODRS
Long	ODSL	ODRL




Multi Abutment Machine Driver

- Dedicated machine driver for Multi Abutment



Abutment Holder 06.2017

- Supplementary instrument for convenient connection of a 2-piece abutment which is difficult to hold with a hand in all oral regions




Abutment Positioning Driver 07.2019

- Used for assembling the abutment in the prosthetic stage after placing a fixture
 - ※ For transfer abutment only
- Function to help convenient and stable mounting and tightening of the abutment kept being pushed away by gingiva
- Used according to the H and G/H lengths of the abutment to be removed as shown below

(Unit : Won)

Range of Use	Short					Long				
	=<9					=>10				
H + G/H	5	6	7	8	9	10	11	12	13	14



Multi Abutment Outer Driver

- Dedicated Torque Driver for Multi Abutment




Locator ☒ Torque Driver

- Dedicated Torque Driver for Locator Abutment



Osstem Torque Driver



- Dedicated driver for osstem torque, which may not be compatible with a general hand piece
- Used after matching the triangle on the outside of the driver with the groove or side of the abutment
- Solid, excellent solid driver only compatible with Ø4.8
- 1.2 hex type L : 5



L \ Type	1.2 Hex	Rigid 4.0	Rigid 4.5	Rigid 5.0	Rigid 6.0	Solid	Excellent Solid
Short (10)	OTH12S	OTR40S	OTR45S	OTR50S	OTR60S	OTS48S	OTE48S
Long (15)	-	OTR40L	OTR45L	OTR50L	OTR60L	OTS48L	OTE48L

Path Probe (TS)

- Checking the path and measuring the gingival height after TS Fixture placement
- C = Connection




Path Probe (KS) 11.2019

- Checking the path and measuring the gingival height after KS Fixture placement
- C = Connection



Torque Connector

- Connector for connecting the torque square driver with a two-way Torque Wrench




Reamer Bite

- Cutting edge to remove lip from the inside of the casted body after casting plastic coping





Machine Driver Connector

- Connector for connecting the machine driver with a two-way Torque Wrench





Reamer Tip (Rigid Abutment)

- Guide part inserted into the casted body for removing lip from the inside after casting plastic coping (for Rigid Abutment)





Driver Handle

- Used by connecting to the Torque Driver



Reamer Tip (Solid, Excellent Solid Abutment)

- Guide part inserted into the casted body for removing lip from the inside after casting plastic coping
- For Solid Ø6.0 and excellent solid Ø4.8
- P= Platform




Finishing Reamer Set

- Used for removing lip from the inside of the casted body after casting plastic coping



Reamer user guide

1. Connected to the casted burn-out cylinder by selecting the reamer tip of the same size as the abutment
2. Rotating the reamer bite with constant force by holding the casted body
3. Reaming until no cutting occurs





For **TSII / III** **SSII / III** **USII / III** **KSIII**

Top panel components

- Hydraulic Membrane Lifter Tube SNMT
- Bone Carrier Head SNBCH30
- Bone Carrier SNBCS35
- Bone Condenser SNBC1120


Lower panel components

- CAS Drill (Ø2.8) SNDR2813TS SNDR2813TL
- CAS Drill (Ø3.1) SNDR3113TS SNDR3113TL
- CAS Drill (Ø3.3) SNDR3313TS SNDR3313TL
- CAS Drill (Ø3.6) SNDR3613TS SNDR3613TL
- CAS Drill (Ø3.8) SNDR3813TS SNDR3813TL
- CAS Drill (Ø4.1) SNDR4113TS SNDR4113TL
- Stopper (7) SNST7 Yellow
- Stopper (2,8) SNST2 SNST8 Purple
- Guide Drill SNGD2027TL
- Twist Drill SNTD2213TL
- Hydraulic Membrane Lifter OCHML
- Stopper (6,12) SNST6 SNST12 Blue
- Stopper (5,11) SNST5 SNST11 Purple
- Stopper (3,9) SNST3 SNST9 Blue
- Stopper (4,10) SNST4 SNST10 Yellow

Crestal Approach – Sinus KIT HÖSSNER

CAS Drill

- Safe lifting of the membrane while forming conical bone for maxillary sinus floor procedure
- Excellent bone removal at low-high speed, and autogenous bone collection at low speed
- Stopper assembled for safe lifting
- Final drill diameter selected based on the bone quality regardless of Straight or Tapered Fixture type
- Recommended speed : 400~800rpm (400rpm for first use)




L \ D	Ø2.8	Ø3.1	Ø3.3	Ø3.6	Ø3.8	Ø4.1
Short	SNDR2813TS	SNDR3113TS	SNDR3313TS	SNDR3613TS	SNDR3813TS	SNDR4113TS
Long	SNDR2813TL	SNDR3113TL	SNDR3313TL	SNDR3613TL	SNDR3813TL	SNDR4113TL

Guide Drill

- Drill to mark the fixture placement position
- Used for removing side walls in a fresh extraction socket with side edges
- Marking line at 2mm from the tip


D	Ø2.0 / 2.7
	SNGD2027TL



Twist Drill (Ø2.2)


- Drilling 1mm under the remaining bone recommended
- Stopper assembled for safe lifting
- End line tip : 0.6mm

D	Ø2.2
	SNTD2213TL



Hydraulic Membrane Lifter Set

- Hydraulic lifting instrument for maxillary sinus membrane
- Winged design with optimized sealing for flapless procedure



Stopper

- Number on the stopper indicating the protruding length of the tip when assembled to a drill or instrument
- Color coded by length
- Drill and stopper use cycle : 50 times




Bone Carrier

- Used for filling the inside of the sinus with bone
- Mounting the head by fastening the back of the body
- Replaceable head (SNBCH30 or SNBCH35) for use



Bone Carrier Head

- Used for filling the inside of the sinus with bone
- SNBCH30 : used after drilling with CAS Drill Ø3.1/3.3
- SNBCH35 : used after drilling with CAS Drill Ø3.6/3.8/4.1 drilling
- Used repeatedly by filling the back of the marking line of the head and taking little by little with a bone condenser to completely fill the inside of the sinus




Bone Condenser

- Instrument to push in the bone material into the sinus




Hydraulic Membrane Lifter Tube

- Connected to the hydraulic membrane lifter



Membrane Lifter 01.2016

- Safe lifting of the membrane due to the round shape with no cutting edge
- lifter selected according to the CAS-Drill diameter as membrane lifting is performed after using the CAS-Drill (head diameter is CAS Drill diameter - 0.2mm)
- CAS Stopper assembled and used for adjusting the depth
- Recommended speed : 400~800rpm (400rpm for first use)
- Be sure to use a drill with irrigation



L \ D	Ø 2.6	Ø 2.9	Ø 3.1	Ø 3.4	Ø 3.6	Ø 3.9
Short	SNML2813TS	SNML3113TS	SNML3313TS	SNML3613TS	SNML3813TS	SNML4113TS
Long	SNML2813TL	SNML3113TL	SNML3313TL	SNML3613TL	SNML3813TL	SNML4113TL



Depth Gauge

- Checking internal lifting of the sinus and measuring the remaining bone depth



Bone Spreader

- Instrument for spreading the filled bone using the engine
- Assembled with a stopper for use
- Recommended speed : ≤ 30rpm (low speed mode)



Y-Connector

- Y-shaped connector for hydraulic lifting of 2 drilling holes at the same time




LAS KIT (HLRSNK) 10.2018




- Lateral Approach - Sinus KIT (LAS KIT) : KIT optimized for lateral approach in maxillary sinus floor procedure
- Including dome drill and core drill for safe formation of a lateral window; and Ø5.5/7.0 diameters available according to the size of the window
- The depth can be adjusted by installing a stopper on the LAS Drill, and the window can be safely formed without perforating the membrane

LAS Full KIT (HLRSNKP) 07.2018

- KIT with 6 additional sinus lift instruments to LAS KIT




LAS KIT Plus Lower Plate



Dome Drill 04.2012

- Forming a window while collecting bone
- Enhanced cutting force with macro and micro cutting edges in combination
- Depth adjusted by assembling with a stopper
- Recommended speed : 1,200~1,500rpm
- ※ Over drilling may result in damage to the membrane

L \ D	Ø5.5	Ø7.0	Wide Ø7.0
25	LSDR554TD	LSDR74TD	LSDR74WTD



Bone Separator 07.2013

- Removing the bone lid from the inside of the core drill



Core Drill 04.2012

- Forming a window while forming the bone lid
- Excellent cutting force and membrane stability due to CAS Drill design concept
- Recommended speed : 1,200~1,500rpm
- ※ Over drilling may result in damage to the membrane

L \ D	Ø5.5	Ø7.0
25	LSDR554TC	LSDR74TC




Stopper 05.2012

- Number on the stopper indicating the protruding length of the tip when assembled to a drill or instrument
- Color coded by length
- Drill and stopper use cycle : 50 times

L	0.5	1.0	1.5	2.0	2.5	3.0
	LSNST0.5	LSNST1.0	LSNST1.5	LSNST2.0	LSNST2.5	LSNST3.0
Color	Yellow	Purple	Blue	Yellow	Purple	Blue

Side Wall Drill 06.2012

- Expanding the window after drilling with a dome drill
- Cutting at 1mm above the lowest part of the drill edge recommended
- Recommended speed : 1,500rpm



Side cutting edge height (mm)	1.0	2.0	3.0	4.0	5.0
CAS KIT Stopper (mm)	8.0	9.0	10	11	12
Side wall drill + CAS KIT Stopper					





※ Depth adjusted by the common CAS KIT Stopper



Lower panel components


For TSII / III SSII / III USII / III KSIII



Crest Remover

- Marking the fixture placement position after removing the narrow ridge horizontally
- Recommended speed
 - Angled type : 1,200~1,500rpm
 - Straight type : 15,000~30,000rpm


L \ D	Ø5.0	Ø7.0
29	CERM50A	CERM70A
45	CERM50S	-



Twist Drill

- Marking the fixture placement position
- Depth adjusted by assembling a stopper according to the fixture length
- Recommended speed : 1,200~1,500rpm


L \ TL \ D	Ø1.8
8.5 \ 33	2D1808LC01
10 \ 34.5	2D1810LC01
11 \ 36	2D1811LC01



Saw 06.2018

- Incision of the narrowed ridge
- After vertical incision, incision in the mesial → distal directions
- Recommended speed : 1,200~1,500rpm
- Recommended use cycle : 10 times
- T = Thickness


T \	Ø7.0	Ø10.0	Ø13.0
0.3	HSAW070	HSAW100	HSAW130



Expansion Drill

12.2016

- Expansion of the ridge after incision
- Used in sequence according to the fixture diameter
F4.0 : I → II → III / F4.5 : I → II → III → IV
- Recommended speed : 25~35rpm



L \ Type	I	II	III	IV
D1 / D2	Ø1.6 / 2.8	Ø2.2 / 3.6	Ø2.8 / 4.4	Ø3.2 / 4.7
8.5	EXP 162808	EXP 223608	EXP 284408	EXP 324708
10	EXP 162810	EXP 223610	EXP 284410	EXP 324710
11.5	EXP 162811	EXP 223611	EXP 284411	EXP 324711

Mount Extension


- Used to exerting torque in manual mode in the process to place or remove an expansion drill into alveolar bone



EXP Mount Driver

- Used to exerting torque for engine in the process to place or remove an expansion drill into alveolar bone

L	
Short (L)	AESMDS
Long (L)	AESMDL



Saw Protector

05.2015

- Safe approach for sawing with a semi-circular saw cover
 - Excellent treatment visibility by forming a window
 - Flexible procedure with a 360° rotating saw
 - Contra angle type (removable saw cover) : KaVo(CL 3-09, S201L), W&H(WS-75)
 - Straight type (built-in saw cover) : KaVo(CL10) ※ Dedicated saw used
- ※ Cover and body of the contra angle type sold separately



Type \ D		Ø7.0	Ø10.0	Ø13.0	Ø15.0	Full Set
Kavo	Contra Angled	Cover Set	SP07AC SP07A	SP10AC SP10A	SP13AC SP13A	- - SP071013A
		Saw Set	- -	SAW10S SP10S	SAW13S SP13S	- SP101315S
	Straight					
W&H	Contra Angled	Cover Set	SP07ACW SP07AW	SP10ACW SP10AW	SP13ACW SP13AW	- - SP071013W

Torque Wrench


- Used for exerting torque on an expansion drill



Depth Gauge

- Instrument to release excessive torque by rotating the hex of the expansion drill with an open wrench when the hand piece does not move with the expansion drill stuck in alveolar bone





Metal Probe

- Instrument to measure the depth of periodontal disease
- Measuring periodontal pockets and identifying the shape of the periodontal pockets such as depth/size
- Marking line for probing in 1 mm increments



HICPM

Plastic Probe

- Instrument to measure the depth of infection or periodontal disease around the implant
- Scratching of implant prevented by using plastic material
- Flexible probe suitable for the curved form of alveolar bone
- Autoclave can be used
- Marking line for probing in 1 mm increments



HICPP

Curette

- Instrument for removing subgingival sediment that is firmly attached to the granulation tissue of a specific area
- Gracey curette
- 01-02 : For removal of granulation tissue from anterior region
- 11-12 : For removal of granulation tissue from the mesial surface in anterior region
- 13-14 : For removal of granulation tissue from the distal surface in anterior region



Type	01-02	11-12	13-14
	HICC0102	HICC1112	HICC1314

Protect Screw

- Preventing infiltration of foreign substances into the internal connection of the fixture using SmartBrush 2
- Tightened with 1.2 hex driver at 5Ncm



Type	Mini	Regular
	OPIB2PM	OPIB2PR

SmartBrush 1

- Used when cleaning Peri-implantitis
- Used after connecting the Protect Screw to the fixture after removing the patients prosthesis or abutment
- Recommended speed : 1,200~1,500 rpm
- Recommended use cycle : About 1 minute per thread
 - ※ Do not use for longer than 4 minutes
- Be sure to polish with saline irrigation and suction
 - ※ Disposable, Do not reuse (Be sure to discard after use)



L	
Short	OPIB1S
Long	OPIB1L

SmartBrush 2 11.2017

- Used for Peri-implantitis cleaning
- Used after connecting the Protect Screw to the fixture after removing the patients prosthesis or abutment
- Be sure to polish with saline irrigation and suction.
- Recommended speed : 1,200~1,500rpm
- Recommended use cycle : 1~2 minutes
 - ※ Excessive use for longer than 3 minutes may result in fracture or bending of the product.
 - ※ Disposable, Do not reuse (Be sure to discard after use)



L \ D	F3.0 / F3.5	F4.0 / F4.5	F5.0 / F5.5	F6.0	F7.0
Short	OPIB23S	OPIB24S	OPIB25S	OPIB26S	OPIB27S
Long	OPIB23L	OPIB24L	OPIB25L	OPIB26L	OPIB27L

Metal Scaler

- Used for removing scale or foreign substance from the surface of the fixture by connecting to an ultrasonic scaler
- Used secondarily after using SmartBrush 1 or SmartBrush 2
- Bendable tip of the product for easy access
- EMS, KaVo and SATELEC types

Type	EMS	KaVo	SATELEC
	HICSME	HICSMK	HICSMS



Plastic Scaler Connector

- Used by assembling to a plastic scaler tip
- Do not use for removing foreign substances from the fixture surface
- EMS, KaVo and SATELEC types
- A = Angle

A \ Type	EMS	KaVo	SATELEC
125°	HICSCE	HICSCK	HICSCS
100°	HICSCE2	HICSCK2	HICSCS2



Plastic Scaler Tip

- Used for removing foreign substances from the abutment or crown by connecting to a SmartScaler
 - ※ Do not use on the fixture surface
- Packing unit : 30ea/1set

1.003.8168



• Including the same components as ESR KIT, which can hold the components provided by other companies

For

Nobel Biocare

Active/Replace /

Straumann

Bone Level /

Astra

Osseo Speed TX



3i Full OSSEOTITE Tapered Certain /

Zimmer

Tapered /

Biohorizons

Internal



ESR Full KIT Surgical Instruments







Not included in the KIT

Guide								
Nobel	Active	Replace	3i	Tapered Certain		Straumann	Bone Level	Roxolid SLActie
	OGNA01L	OGNR02L		OGIF01L	OGS801L		OGSTR5	
	OGNA02L	OGNR03L		OGIF02L	OGS802L		OGSTRL	
		OGNR04L						
Astra	Osseo Speed TX		Biohorizons	Internal	External	Zimmer	Tapered	
	OGA001L			OGZB01L	OGBES		OGZB01L	
	OGA002L			OGZB02L	OGBEL		OGZB02L	
SR Drill		SR Tip	Screw Holder			Re-Tap		
OSRD09L		ORT18L	OSHR18L			ORTR18L		





Guide

- Centering and shaking prevention of SR Drill, SR Tip, etc. by connecting and fixing to the fixture
- Use according to fixture type and diameter
(Internal, submerged type products of 6 overseas companies)
- Short or Long types selected according to the intermaxillary distance
- ■ Used in common
- C = Connection / F = Fixture



Osstem

C \ Type	TS		SS		US		KS	
	Short	Long	Short	Long	Short	Long	Short	Long
Mini Regular Wide							-	-
	OGTRS	OGTRL	OGSRS	OGSRL	OGURS	OGURL	OKGRS	OKGRL
	-	-	OGSRS	OGSRL	OGUWS	OGUWL	-	-
	-	-	-	-	-	-	-	-





Nobel Biocare

F \ Type	Active		Replace	
	Short	Long	Short	Long
Ø3.5 Ø4.3 Ø5.0 Ø6.0				
	OGNA01S	OGNA01L	-	-
	OGNA02S	OGNA02L	OGNR02S	OGNR02L
	OGNA02S	OGNA02L	OGNR03S	OGNR03L
	-	-	OGNR04S	OGNR04L



Nobel Biocare

F \ Type	MkIII	
	Short	Long
Ø3.3 Ø3.75 Ø4.0 Ø5.0		
	OGUMS	OGUML
	OGURS	OGURL
	OGURS	OGURL
	OGUWS	OGUWL



Straumann

F \ Type	Bone Level		F \ Type	Roxolid SLActive	
	Short	Long		Short	Long
NC (3.3) RC (4.1) RC (4.8)			RN (3.3) RN (4.1) RN (4.8) WN (4.8)		
	OGSB01S	OGSB01L		OGSTRS	OGSTRL
	OGSB02S	OGSB02L		OGSTRS	OGSTRL
	OGSB02S	OGSB02L		OGSTRS	OGSTRL


Astra

F \ Type	Osseo Speed TX	
	Short	Long
Small (3.5 S) Small (4.0 S) Large (4.5) Large (5.0) Large (5.0 S)		
	OGAO01S	OGAO01L
	OGAO02S	OGAO02L
	OGAO02S	OGAO02L
	OGAO02S	OGAO02L
	OGAO02S	OGAO02L



3i

F \ Type		Full Osseotite Tapered Certain		F \ Type		Full Osseotite Tapered	
		Short	Long			Short	Long
							
3.25		OGIF01S	OGIF01L	Ø4.0		OGURS	OGURL
4.0		OGIF02S	OGIF02L	Ø5.0		OGURS	OGURL
5.0		OGIF02S	OGIF02L	Ø6.0		OGURS	OGURL
6.0		OGIF02S	OGIF02L				

Zimmer

F \ Type		Tapered	
		Short	Long
			
Green (3.7)		OGZB01S	OGZB01L
Green (4.1)		OGZB01S	OGZB01L
Green (4.7)		OGZB01S	OGZB01L
Green (6.0)		OGZB02S	OGZB02L

Biohorizons

F \ Type		Internal (Tapered Bone Level)		F \ Type		External	
		Short	Long			Short	Long
							
Yellow		OGZB01S	OGZB01L	Ø3.5		OGUMS	OGUML
Green		OGZB01S	OGZB01L	Ø4.0		OGURS	OGURL
Blue		OGZB02S	OGZB02L	Ø5.0		OGBES	OGBEL
				Ø6.0		OGBES	OGBEL

Reverse Drill 06.2017

- Instrument used for removing fractured screws
- Be sure to use with a suitable guide for the fixture
- When the red marking of the reverse driver is shown above the guide assembled to the fixture, use a screw holder to remove the fractured screw
- For hand mode / Rotating direction : Reverse rotation / Use cycle : 10 times
- ※ Do not use more than 10 times. Do not reuse
- F = Fixture

L \ Type	M1.6	M1.8	M2.0
Short	-	ARVDRSC	ARVDRSC
Long	ARVDMLC	ARVDRLC	ARVDRLC



Screw Removal Drill (SR Drill) 12.2014

- Used for removal to form a hole in fractured screws
- Be sure to assemble the guide and remove the cut chips by suction with irrigation into the window
- Short and Long types according to the intermaxillary distance
- Drilling until the red line around the handle is not visible
- Recommended speed : Reverse rotation of 1,200~1,500rpm / Use cycle : 5 times
- ※ Be sure to use with a guide assembled. / Do not exert excessive vertical force. / Do not soak in hydrogen peroxide.
- ※ Disposable, Do not reuse
- Short : Sold separately

L \ Type	M1.6	M1.8	M2.0
Short	OSRD08S	OSRD09S	OSRD10S
Long	OSRD08L	OSRD09L	OSRD10L



Torque Driver Handle

- Used by rotating with a hand after assembling with products such as SR tip, AR tip, and screw holder



Reverse Driver

10.2010

- Instrument used for removing fractured screws
- Be sure to use with a suitable guide for the fixture
- When the red marking of the reverse driver is shown above the guide assembled to the fixture, use a screw holder to remove the fractured screw
- For hand mode / rotating direction : reverse rotation / use cycle : 10 times
- ※ Do not use more than 10 times
- F = Fixture

L \ F	Mini	Regular /Wide
Short	-	ORVDRS
Long	ORVDML	ORVDRL



Re-tap

- Instrument to restore the thread to the initial state when the screws cannot be fastened due to damage to the internal thread of the fixture
- Thread formed in hand mode with a Torque Wrench or ratchet wrench

Type	M1.6	M1.8	M2.0
	ORTM	ORTR18	ORTR



Screw Removal Tip (SR Tip)

- Removing fractured screws by rotating the screw removal tip in the hole in the fractured surface of the screws formed by using the screw removal drill(SR Drill)
- Rotating direction : Reverse rotation
- ※ Disposable, Do not reuse

L \ Type	M1.6	M1.8	M2.0
Short	ORT16S	ORT18S	ORTS
Long	ORT16L	ORT18L	ORTL



ESR Handle

03.2013

- Instrument to fix the guide to the fixture

Type	OARH
------	------



Screw Holder

- Removing partially protruding fractured screws by assembling with a screw holder
- Color coded for easy type indication
- Rotating direction : Reverse rotation

Type	M1.6	M1.8	M2.0
	OSHM	OSHR18	OSHR



Abutment Removal Tip (AR Tip)

07.2017

- Used for partial fractured abutment, mount remaining and stuck in the fixture
- After assembling it to the fractured abutment hole and fixing in place, remove by shaking with a forcep, etc.
- Mini : removing screws with a slipped hex
 - After assembling it to the slipped hex, rotate in the reverse direction to connect to the screw for removal

L \ Type	Mini	Regular
Short	OARTMS	OARTRS
Long	OARTML	OARTRL
Ex.Long	OARTMEL	OARTREL



Slot Driver

10.2010

- Instrument to use by forming a slot with Ø0.8 bur, when force cannot be exerted using a driver due to the damaged hex of Healing Abutment, Cover Screw, or Abutment Screw

OTSD07



Transfer Abutment Separate Tool

01.2009

- Used to release the jamming caused by Non-hex Transfer Abutment stuck due to the contact of the fixture and the morse taper
- Common use, by using the body end for mini, placing regular into the 2nd groove
- It is easy to remove if the body and abutment are integrated by rotating the driver forward after removing the abutment screw and placing a separate tool body into the inner hole of the abutment

If it is difficult to separate, use after connecting a ratchet wrench to the driver

Driver	Body	Set
TASD	TASB	TAST



Driver



Body



EFR KIT

Easy Fixture Removal KIT (OSFRK)

RENEWAL 2020



Top panel components

Fixture Wrench
FRDFE



For **TSII / III** **SSII / III** **USII / III** **KSIII** **Ultra-wide**

EFR Full KIT

Easy Fixture Removal Full KIT (OSFRFK)

01.2018

- Including the same components as EFR KIT, which can hold the components provided by other companies

Lower panel components

Fixture Wrench
FRDFE





Torque Wrench
TW400B



For

Nobel Biocare	Active/Replace	/	Straumann	Bone Level	/	Astra	Osseo Speed IX
3i	Full OSSEOTITE Tapered Certain	/	Zimmer	Tapered	/	Biohorizons	Internal



EFR Full KIT Surgical Instruments

Not included in the KIT

Remover Screw							
Nobel	Active			Replace			
	Normal	Fracture		Normal	Fracture		
	FRSMNA35 FRSR40 FRSW50	FRSMNA35F FRSR40F FRSW50F		FRSMNR35 FRSR40 FRSW50	FRSMNR35F FRSR40F FRSW50F		
Straumann	Bone Level		3i	Full Osseotite Tapered Certain		Biohorizons	Internal
	Normal	Fracture		Normal	Fracture		
	FRSM33 FRSR541 FRSW548	FRSM33F FRSR541F FRSW548F		FRSMI325 FRSRI40 FRSWI50	FRSMI325F FRSRI40F FRSWI50F		Normal FRSRZ41 FRSWZ47 FRSWZ60
Zimmer	Tapered		Astra	Osseo Speed TX		Remover Body	
	Normal	Fracture		Normal	Fracture		
	FRSMZ37 FRSRZ41 FRSWZ47 FRSWZ60	FRSMZ37F FRSRZ41F FRSWZ47F FRSWZ47F		FRSMNA35 FRSRA40 FRSR40 FRSW50	FRSMNA35F FRSRA40F FRSR40F FRSW50F		FRBW57S FRBW57L FRBUW60S FRBUW60L

Remover Screw

- Acting as a support structure for reverse rotation of the remover body after connected and fixed to the fixture
- Used according to the type and diameter of the fixture to remove (Internal/submerged type products of 6 overseas companies, normal/fracture)
- Fracture used for removing fixtures with the hex entirely fractured
- Compatible with products of 6 overseas companies
- Recommended tightening torque : Regular/Wide 80Ncm, Mini 60Ncm
- T = Type ※ Disposable, Do not reuse



Osstem

T \	Mode	Mini Ø3.5 /-	Regular Ø4.0~4.5 / P4.8	Wide Ø5.0 / P6.0
TS/SS	Normal	FRSM35	FRSR40	FRSW50
	Fracture	FRSM35F	FRSR40F	FRSW50F
US		FRSM35US	FRSR40US	FRSW50US
KS	Normal	KSFRSM35	KSFRSR40	KSFRSW50
	Fracture	KSFRSM35F	KSFRSR40F	KSFRSW50F

Nobel Biocare

T \	Mode	Mini Ø3.5	Regular Ø4.3	Wide Ø5.0/6.0
Active	Normal	FRSMNA35	FRSR40	FRSW50
	Fracture	FRSMNA35F	FRSR40F	FRSW50F
Replace	Normal	FRSMNR35	FRSR40	FRSW50
	Fracture	FRSMNR35F	FRSR40F	FRSW50F

Straumann

T \	Mode	Mini Ø3.3	Regular Ø4.1	Wide Ø4.8
Bone Level	Normal	FRSMS33	FRSR541	FRSW548
	Fracture	FRSMS33F	FRSR541F	FRSW548F

Astra

T \	Mode	Mini Ø3.5	Regular Ø4.0	Regular Ø4.5	Wide Ø5.0
Osseo Speed TX	Normal	FRSMNA35	FRSRA40	FRSR40	FRSW50
	Fracture	FRSMNA35F	FRSRA40F	FRSR40F	FRSW50F

3i

T \	Mode	Mini Ø3.25	Regular Ø4.0	Wide Ø5.0/6.0
Full Osseotite Tapered Certain	Normal	FRSMI325	FRSRI40	FRSWI50
	Fracture	FRSMI325F	FRSRI40F	FRSWI50F

Zimmer

T \	Mode	Mini Ø3.7	Regular Ø4.1	Wide Ø4.7	Ultra-wide Ø6.0
Tapered	Normal	FRSMZ37	FRSRZ41	FRSWZ47	FRSWZ60
	Fracture	FRSMZ37F	FRSRZ41F	FRSWZ47F	FRSWZ47F

Biohorizons

T \	Mode	Mini Ø3.8	Regular Ø4.6	Wide Ø5.8
Internal	Normal	FRSRZ41	FRSWZ47	FRSWZ60
	Fracture	FRSRZ41F	FRSWB46F	FRSWB46F

Screw Driver

- Driver to connect and fix the remover screw to the fixture
- Recommended remover screw tightening torque : Regular/Wide 80Ncm, Mini 60Ncm
- F = Fixture

F	Mini	Regular	Wide
	FRSDM23	FRSDR25	FRSDW30



Torque Wrench

- Used to remove the fixture with the remover body after tightening with screw driver
- Torque applied up to 400Ncm (80/100/200/300/400Ncm scale display)
- Torque applied by aligning the center of the bar with the torque value to be applied by pulling the bar
- Washed and sterilized after use for storing

	TW400B
--	--------



Remover Body

- Instrument to exert torque in the fixture loosening direction by connecting to a remover screw
- Used according to the diameter of the fixture to remove
- ※ Disposable, Do not reuse
- F = Fixture

F	Mini	Regular	Only for osstem Wide	Only for overseas companies Wide	Ultra-wide
Short	FRBM35S	FRBR40S	FRBW50S	FRBW57S	FRBUW60S
Long	FRBM35L	FRBR40L	FRBW50L	FRBW57L	FRBUW60L



Fixture Wrench

- Wrench to remove fixture from the remover body

	FRDFE
--	-------



Torque Extension

- Screw driver and remover body length extension (up to 10mm)

	OTE
--	-----



Dr. Cho's Instrument KIT

(DCHOKIT)

11.2017

- Optimal implant surgery KIT based on years of clinical know-how
- Consisted of 10 types of instruments (1ea each)

Periosteal Elevator (24G)

- Lifting mucosal periosteum after gingival tissue incision
- W : 4.2/4.0mm

EP24G-W

Minesota Retractor

- Securing a clear view by pulling the mouth, cheeks, etc.

RTCRM-W

Extension Hose Adapter

- Adapter for chair suction connection

SNKHA-W

Periosteal Elevator (Selden)

- Lifting mucosal periosteum after gingival tissue incision
- W : 10/13mm

EP23-W

Needle Holder (Crile-Wood, TC)

- Straight
- Tungsten carbide-treated beak
- L : 150mm (±5)

NHC150TC-W

Periodontal Chisel

- Bone removal and formation
- W : 5.0mm

CHCO2-W

Tissue Forcep (ADSON)

- Used to hold soft tissue
- No projections inside the beak
- L : 120mm (±5)

PT41-W

Extension Hose

- Extension hose for chair suction connection
- Autoclave can be used
- Transparent silicone material

SNKHS-1-W

Dr.Cho's Instrument Pouch

- Used for storing and sterilization of instruments
- L : 550 X 400mm

WPB-W

Titanium Suction Tip

- D(Inner Diameter) : 3.0mm

SN3TI-W



Osstem Basic Instrument KIT

(OBKIT)

11.2017

- Universally used implant surgery KIT
- Consisted of 25 types of instruments (1ea each)

Periosteal Elevator (24G)

- Lifting mucosal periosteum after gingival tissue incision
- W : 4.2/4.0mm

EP24G-W

Scalpel Handle (Flat Type)

SHF-W

Needle Holder

- Mayo-Hegar
- Tungsten carbide-treated beak
- L : 160mm (±5)

NH160TC-W

Pouch

- Used for storing and sterilization of instruments
- L : 470 X 400mm

WPA-W

Mirror

MHC-DMCS4-W

Chisel

- Bone removal and formation
- Oschenbien & fedi (curved)
- W : 5.0mm

CHCO2-W

Titanium Suction Tip

- D(Inner Diameter) : 3.0mm

SN3TI-W

Hemostats

- Mosquito (curved)
- L : 130mm (±5)

HTM130C-W

Tissue Forcep ADSON

- No projections inside the beak
- L : 120mm (±5)

PT42-W



Osstem Basic Instrument KIT

(OBKIT) 11.2017

Periosteal Elevator (Selden)

- In case of Genival tissue flap, retract and fix
- W : 10/13mm

EP23-W

Caliper

- Castroviejo

LPC90-W

Bone Well

BWSUS1-W

Scissor (LaGrange)

- Compound (curved)
- L : 115mm (±5)

SCLC115-W

Surgical Curettes (Gracey)

CGR11-12-W

Surgical Curettes (Surgical Curettes, CM11)

URCM11-W

Minesota Retractor

RTCRM-W

Towel Clamp

- Towel Clamp, Backhaus
- L : 135mm (±5)

CPTC135-W

Periosteal Elevator (MOLT9)

- Lifting mucosal periosteum after gingival tissue incision
- W : 10/13mm

EP9-W

Scalpel Handie (Straight Type)

SHS-W

Scissor (Tissue Scissor)

- Straight
- L : 150mm (±5)

SCTC115-W

Mallet

- Autoclave can be used

ML25-W

Tweezer (Wide)

- L : 155mm

PCW150-W

Bone Rongeurs

- Friedman
- L : 140mm (±5)

RNGF140-W

Surgical Curettes (Surgical Curettes, CM10)

URCM10-W


Periosteal elevator (Prichard)

- W : 11/4.9mm

EP9R3-W

Custom KIT (OCTK) 01.2009

- KIT used to disinfect some of the surgical instruments or to store new spare tools
- Additional 3 types of rubber (large, medium, small) which can be used according to user preference
- Sterilizable material (132°C, 15 minutes)




Healing Case (OHAC) 02.2018

- Case for temporary storage and cleaning of Healing Abutment during the prosthesis procedure
 - Upper prosthesis for additional mounting : Transfer / Temporary / Angled / Cover Screw / Pick-up & Transfer Impression coping / OB Anchor / temporary crown (Only the Healing Abutment can be combined with the top plate.)
- Like the tooth arrangement, a total of 28 cells are composed of 7 cells each in the upper / lower and left / right sections
- Sterilizable material (132°C, 15 minutes), sterilization required for reusing the case
- ※ This product is not a case for reuse of Healing Abutment




Osteo KIT (OSTK) 01.2009

- KIT used for maxillary sinus floor elevation to vertically increase the amount of alveolar bone available in the maxillary anterior region
- Expanding osteotome : KIT used to increase the initial fixation stability of the implant by densifying the trabeculae of bone while preserving the bone instead of removing it from low quality bones
- Stopper for adjusting the depth of procedure



Osteotome KIT (AOST) 09.2011


- KIT used for maxillary sinus floor elevation to vertically increase the amount of alveolar bone available in the maxillary anterior region
- Concave type only
- Stopper for adjusting the depth of procedure



Osteotome Stopper 05.2018


- Stopper for procedure depth adjustment, sold separately

D	Ø 2.0	Ø 2.5	Ø 3.0	Ø 3.5	Ø 4.0
	OST20SH	OST25SH	OST30SH	OST35SH	OST40SH





Sinus KIT (ASLK) 01.2009

- KIT containing various tools for maxillary sinus floor elevation (sinus procedure)
- Lateral approach instrument for sinus
- Components (5 types)
 - Freer elevator : OFE
 - Bone Graft Carrier : OBGC
 - Membrane Separator (Circle type) : OMSC
 - Sinus Currette-Short : OSCS
 - Sinus Currette-Long : OSCL




Bone Spreader KIT (OBSOK) 01.2009

- KIT used for expanding narrowed alveolar ridge
- Offset type convenient for surgery
- Components (4 types)
 - OBSO22F, OBSO28F, OBSO35F, OBSO35R



- Use for alveolar bone expansion
- Offset type for easy operation
- Depth marking corresponding to the implant length



		(Unit : mm)			
Code	Tip length	7	8.5	10	11.5
	Spec.				
OBSO22F	Thickness	1.15	1.3	1.45	1.6
	Width	2.1	2.2	2.2	2.2
OBSO28F	Thickness	1.15	1.3	1.45	1.6
	Width	2.65	2.8	2.8	2.8
OBSO35F	Thickness	1.3	1.45	1.6	1.8
	Width	3.3	3.5	3.5	3.5
OBSO35R (round type)	Thickness	1.85	2.1	2.3	2.55
	Width	3.3	3.5	3.5	3.5




Ridge Split KIT

Straight (ORSSK)

01.2009

Straight

- Chisel : Used for expanding narrowed alveolar ridge
- Blade Holder : Malleting enabled by tightening a #15 blade when it is difficult to make a bone incision using bur due to low bone quality
- Components
 - Ridge Split Chisel : ORSS15, ORSS20, ORSS25, ORSS30
 - Blade Holder : ORSBH



Code	Spec.	Tip length			
		7	8.5	10	11.5
ORSS15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSS20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSS25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSS30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4




Ridge Split KIT

Offset (ORSOK)

01.2009

Offset

- Chisel : Used for expanding narrowed alveolar ridge
- Blade Holder : Malleting enabled by tightening a #15 blade when it is difficult to make a bone incision using bur due to low bone quality
- Components
 - Ridge Split Chisel : ORSO15, ORSO20, ORSO25, ORSO30
 - Blade Holder : ORSBH



Code	Spec.	Tip length			
		7	8.5	10	11.5
ORSO15	Thickness	1.1	1.27	1.5	1.5
	Width	4	4	4	4
ORSO20	Thickness	1.45	1.7	2.0	2.0
	Width	4	4	4	4
ORSO25	Thickness	1.8	2.15	2.5	2.5
	Width	4	4	4	4
ORSO30	Thickness	2.15	2.5	3.0	3.0
	Width	4	4	4	4

Osstem Implant System product description

Osstem Implant offers a variety of dental bone graft materials, as well as fixtures made of medical grade titanium. Osstem Implant's abutments, prosthetic materials and surgical tools are only compatible with Osstem Fixtures. If used with products of other manufacturers, it may cause problems including loosening and fractures due to incomplete tightening and compatibility. For more details about any individual product, please refer to the user manual, catalog or visit our company website (www.osstem.com). Please check the product labels for product codes, specifications, date of manufacture and expiration date.

Sterilization

Fixtures, cover screws and Healing Abutments are pre-cleaned and sterilized by gamma rays. These products are sterile, disposable medical instruments and must be handled in a sterile field using sterilized tools to prevent contamination and infection of the product or treatment area. If the package has been opened, damaged or has expired, the product must be discarded due to the risk of contamination, infection and bone graft failure. If re-sterilized or re-used, the product may result in infection, osseointegration failure, and damage to implants due to reduced precision.

Storage conditions

Store in a dry place at room temperature (1~30°C). Keep away from direct sunlight.

General precautions

Warning! Implant surgical techniques involves professional and complex processes. To perform dental implant surgery, relevant formal training and education is required. If the patient has bone disease (osteoporosis, osteomalacia) or metabolic bone disorders, special considerations should be given to these conditions prior to surgery.

Precautions

Suitability of bone and proper surgical procedures should be taken into account when determining an implant surgery. Proper implant should be prepared in consideration of anticipated situations and precautions. Excessive occlusal load may cause loosening or fracture of an implant. In order to avoid this condition, the implant must be placed in accurate location and direction considering the relationship between the implant and opposing dentition. Visual inspection as well as radiographic examinations are essential to determine basic presurgical information, occlusal conditions and adequacy of the bone. Adequate radiographs, surgical planning and visual inspection of the implant site are required prior to implant surgery.

Procedural precautions

Osstem Implant System is for single- or two-stage procedure. Special attention should be paid to temperature, surgical lesions and removal of the sources of contamination and infection in an attempt to minimize damage to the cell tissue. All drills and taps must be continuously and sufficiently irrigated for cooling. Implant placement should be accomplished at very low speed (25~30 rpm) or manually. Excessive torque (greater than 55Ncm) can have adverse effects such as partial fracture or necrosis of the bone. Placing an implant tilted by 30° or higher is not recommended due to possible fracture or implant. Immediate loading to the fixture right after the surgery should be avoided. The bone quality and initial stability after fixture placement are important elements in determining the appropriate loading time. Mini-diameter implant or implant with diameter of 4.0 or less which integrates with Angled Abutment may be fractured due to limitations of structural rigidity. They are not recommended for use in a posterior area. Ultra-wide Fixtures are intended to be used only in the posterior region and should not be used with Angled Abutments.

If considering the Ultra-wide fixtures, radiographic evaluation should be performed to determine the bone mass and potential anatomical restrictions. Short implants (diameter greater than 5mm, shorter than 7mm) are only used for the posterior region. Clinicians must thoroughly examine the patient for any of the following conditions: 1) Peri-implant bone loss, 2) Changes to implant's response to percussion, 3) Vertical changes in the osseointegrated fixtures determined by X-ray. If a short implant shows loosening or greater than 50% bone loss, the implant should be considered for possible removal. Clinicians should consider a two-stage surgical approach, splinting with other implants and placement of the widest possible diameter fixture. Allow sufficient healing time for osseointegration before prosthesis and avoid immediate loading. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load. Avoid applying HA-coated fixtures to hard bone because damage and cracks might occur in the coated layer. It is recommended that the insertion torque of the implant be less than 35Ncm. The surfaces of CA and SOI have the same physical shape as the SA surface made through blasting and etching treatments. These surfaces are designed to maintain the SA surface chemically-activated by encasing CA in a solution and SOI in a hydrophobic coating after the SA surface treatment to prevent the product from being exposed to air. Thus, CA or SOI products should be placed in the target region at least within 15 minutes after removal from the vial.

Warning

Improper patient selection and treatment planning may result in dental implant failure or loss of bone supporting the implant. Osstem Implant System must not be used for purposes other than intended and must not be altered in any shape or form. Implant loosening, bone loss and chronic infections can result in implant failure.

Indications for use

Osstem Implant System is an artificial dental root that has been designed for use in dental implant treatment for restoring missing teeth. It can be placed via surgical procedures in maxillary or mandibular bone to replace natural dental root. The System is intended for use in fabricating temporary or final appliances in the form of cement-retained, screw-retained, overdenture and fixed-bridge to replace a single tooth or multiple teeth in the maxillary/mandibular region or for partially or fully edentulous patients. Products with diameter of 3.25mm or less must be used exclusively for mandibular anterior teeth in order to prevent fracture due to excessive occlusal load.

Side effects

There are possible side effects after implant surgery (loss of implant stability, damaged prosthesis, etc.). These issues can be caused by the lack of available bone or poor bone quality, infection, patient's poor oral hygiene or non-compliance with post-op procedures, allergic reaction, movement of the implant, degradation of surrounding tissue, or improper placement/arrangement of the implant.

Contraindications

- Contraindications include the following, but are not limited to:
- Patients with hemophilia or issues related to bone or wound treatment
 - Patients with uncontrollable diabetes or patients that smoke or drink excessively
 - Patients with compromised immune systems due to disease or chemo/radiation therapy
 - Patients with oral infection or inflammation (improper oral hygiene or bruxism)
 - Patients with incurable malocclusion/joint disorder and insufficient dental arch space
 - Patients who are not suitable for surgery.



■ Manufacturer : Osstem Implant Co., Ltd.
203, Geoje-daero, Yeonje-gu, Busan, Korea
TEL 82-51-850-2500 FAX 82-51-861-4693



DEUTSCHE OSSTEM GmbH.
Mergenthalerallee 25
65760 Eschborn, Germany
+49-(0)6196-777-550

Storage condition
Dry place at room temperature

Rx only
For USA only : Federal law restricts this device to sale by or on the order of a dentist



OSSTEM[®]
IMPLANT