

Sinus Membrane Elevation without Perforation


CAS KIT

- Prevents perforation of sinus membrane
- Easy and safe membrane elevation with hydraulic pressure
- Applicable to various cases




Prevents Perforation of Sinus Membrane

- Prevents perforation of sinus membrane and gets rid of patient pain by using round shaped drill without malletting





Osteotome KIT





CAS KIT

Easy and Safe Membrane Elevation with Hydraulic Pressure

- Hydraulic pressure lifts membrane widely and high.
- Elevated sinus has excellent blood supply and osteogenesis capability




Round shape
safe drilling




Creates inverted conical
shaped bone lid
Prevents perforation


Applicable in Various Cases



Flat



Septum



Inclined



Nerve