

Dr. Cho's instrument Kit

Order Code : DCHO KIT

Optimal Implant Surgery Instrument Kit with Long-term Clinical Know-how



OSSTEM[®]
IMPLANT

1. Periosteal Elevator

- Comparing to Molt 9, with narrow width, fine and precise surgery possible
- The instrument lifting the mucosal periosteum after cutting the gingival tissue



2. Minesota Retractor

- Improving gripping area for long time use, which reduces hand fatigue (patent)
- The device that can secure the visibility by pulling the mouth, cheeks, etc.



3. Periosteal Elevator Selden

- The device that lifts and fixes the mucosal periosteum after cutting the gingival tissue




4. Needle Holder, Crile-Wood, TC

- Preventing the situation of being caught intermittently by bevel processing in the middle of the needle holder
- Tungsten carbide treatment in the Beak area enables retaining fine sutures, thereby increasing the convenience of the procedure.




5. Ochsenbein Chisel

- Useful when collecting a small amount of autogenous bone
- The blade is toward the front




6. Tissue Forcep, ADSON

- The instrument catching the gingival without hurting the soft tissue due to no protrusion on the inner surface.



7. Titanium Suction Tip

- High durability with titanium tip (inner diameter 3mm)



8. Dr.Cho's Instrument Pouch

- Using for the storage and sterilization of instruments
- Size 470X400mm

9. Extension Hose

- Hose for extension, connecting to unit chair suction
- Possible for autoclave sterilization (silicone material)
- Cleaning in transparent, Easy in dealing with

10. Extension Hose Adapter

- Adaptor to connect to unit chair suction

Osstem Basic Instrument Kit Product Code : OB KIT

The Kit with 25 Kinds of Instruments Composed of Commonly used Surgical Instruments

Periosteal Elevator



Selden



MOLT9



Prichard



P24G

Mirror



Mirror

Scissor



Tissue scissor



Scissor, LaGrange

Curette



Gracey



Surgical Curette, CM10



Surgical Curette, CM11

Hemostat



Mosquito

Needle Holder



Mayo-Hegar

Towel Clamp



Backhaus

Tissue Forcep



ADSON (ToothO)

Tweezer




Wide Tweezer

Scalpel Handle



Scalpel Handle (straight)



Scalpel Handle (Flat)

Caliper



Castroviejo

Retractor



Minesota Retractor

Chisel



Periodontal Chisel

Rongeurs



Bone Rongeurs, Friedman

Mallet



Mallet

Bone Well



Bone Well

Suction Tip



Titanium Suction Tip

Instrument Pouch



NO.	Code	Description
1	EP9-OB	Periosteal Elevator Molt9
2	EPPR3-OB	Periosteal Elevator Prichard
3	EP23-OB	Periosteal Elevator Selden
4	HTM130C-OB	Hemostat 12cm (Curved)
5	PT42-OB	Tissue Plier, ADSON(SAW TYPE)12cm
6	SCTC115-OB	Iris Scissor T/C
7	SCLC115-OB	La grange Scissors 11cm
8	RNGF140-OB	Bone Rongeur 14cm
9	ML25-OB	Mallet ML-01
10	CGR11-12-OB	Gracey Curette
11	RTCRM-OB	Minesota Retractor
12	NH160TC-OB	Needle Holder T/C 16cm
13	SHS-OB	Scalpel Handle (Straight)
14	SHF-OB	Scalpel Handle (Flat)
15	CHCO2-OB	Periodontal Chisel Ochsenbein & Fedi Curved
16	SN3TI-OB	Titanium Suction Tip
17	PCW150-OB	Wide Tweezer
18	DMCS4-OB	Mirror Set
19	LPC90-OB	Castro Gauge (Curved)
20	EP24G-OB	Periosteal Elevator P24G
21	URCM10-OB	Surgical Curette Curved CM10
22	URCM11-OB	Surgical Curette Curved CM11
23	CPTC135-OB	Towel Clamp
24	BWSUS1-OB	Bone Well (Small)
25	WPA-OB	Instrument Pouch