


## One Touch, Automatic Pressure Adjustment

# SlowJec

Minimize patient pain with automatic pressure adjustment function, excellent usability with one-touch operation

### Automatic Pressure Adjustment by Sensor

- Pause when over pressure is detected by measuring the pressure loaded on the needle
- When the pressure decreases, automatically inject the remaining amount



Pause when over pressure is detected

Automatic injection when the pressure is decreased

### Normal Syringe vs SlowJec

	Normal syringe	SlowJec
Injection time	With in 10seconds	<b>1min 4seconds</b>
Pressure	High	<b>Low</b>
Pain	Strong	<b>Very weak</b>

※ Criteria : Inject 100% 1.8ml, in one touch mode

**LCD Screen**  
Check Injection quantity, pressure quantity, injection speed

**Voice Support**  
Beeps to notify

**AS** **Aspiration Function**  
Check the penetration of the vessel before injection

**One Touch Button**

**Compatible with General Types of Needles**  
Wide range of drugs can be injected

**Replaceable Battery**  
Removable

**Various Injection Volume Control and Speed Adjustment**  
5-level injection volume controls

### Easy One Touch Operation

- Press once to complete the injection
- No need of gripping the syringe tight for the whole time



※ Check Injection quantity, pressure quantity, injection speed

### Convenient Pen Type & Cordless

- Easy to use like a pen · Compact size
- Light weight (150g including a battery)