


Clear Image! Light Weight!

N1 Portable X-ray

- Clear image
- Light weight (approximately 1.5kg)
- High efficiency battery

Clear Image

- As clear as standard X-ray (70 kV, 3mA)
- Clarity enhanced by FOCAL SPOT (0.3mm)



[70kV, 3mA]

[70kV, 2mA]

※ image above is subject to change based on sensor.

Operating with one hand

- No strains on female operator's wrist



Intuitive screen and easy to use

- Intuitive LCD screen displaying only pertinent information (such as target location and person, working time)
- Press buttons to prevent device malfunction and breakdown



Charge only once a week! (given 30 images taken per day)

- Can take 300 images with high efficiency battery

